

CLEC1B Antibody

Catalog No: #35685

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	CLEC1B Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Antigen affinity purification.
Applications	IHC
Species Reactivity	Hu
Specificity	The antibody detects endogenous levels of total CLEC1B protein.
Immunogen Type	Recombinant Protein
Immunogen Description	Full length fusion protein
Target Name	CLEC1B
Other Names	CLEC2; CLEC2B
Accession No.	Swiss-Prot#: Q9P126NCBI Gene ID: 51266Gene Accssion: BC029554
Uniprot	Q9P126
GeneID	51266;
Concentration	0.5mg/ml
Formulation	Rabbit IgG in pH7.4 PBS, 0.05% NaN ₃ , 40% Glycerol.
Storage	Store at -20°C

Application Details

Immunohistochemistry: 1:10-1:50

Images



Immunohistochemical analysis of paraffin-embedded Human liver cancer tissue using #35685 at dilution 1/15.

Background

Natural killer (NK) cells express multiple calcium-dependent (C-type) lectin-like receptors, such as CD94 (KLRD1; MIM 602894) and NKG2D (KLRC4; MIM 602893), that interact with major histocompatibility complex class I molecules and either inhibit or activate cytotoxicity and cytokine secretion. CLEC2 is a C-type lectin-like receptor expressed in myeloid cells and NK cells.

Note: This product is for in vitro research use only