

ENPP5 Antibody

Catalog No: #35725

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	ENPP5 Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Antigen affinity purification.
Applications	IHC
Species Reactivity	Hu
Specificity	The antibody detects endogenous levels of total ENPP5 protein.
Immunogen Type	Recombinant Protein
Immunogen Description	Fusion protein corresponding to a region derived from internal residues of human ectonucleotide pyrophosphatase/phosphodiesterase 5 (putative)
Target Name	ENPP5
Other Names	NPP-5; KIAA0879; ENPP5
Accession No.	Swiss-Prot#: Q9UJA9NCBI Gene ID: 59084Gene Accssion: BC027615
Uniprot	Q9UJA9
GeneID	59084;
Concentration	1.9mg/ml
Formulation	Rabbit IgG in pH7.4 PBS, 0.05% NaN ₃ , 40% Glycerol.
Storage	Store at -20°C


Application Details

Immunohistochemistry: 1:50-1:200

Images



Immunohistochemical analysis of paraffin-embedded Human liver cancer tissue using #35725 at dilution 1/40.



Immunohistochemical analysis of paraffin-embedded Human tonsil tissue using #35725 at dilution 1/40.

Background

ENPP5 is a single-pass membrane protein that belongs to the nucleotide pyrophosphatase/phosphodiesterase family. It may play a role in neuronal cell communication, and it lacks nucleotide pyrophosphatase and lysopholipase D activity.

Note: This product is for in vitro research use only