

FOXR2 Antibody

Catalog No: #35740

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	FOXR2 Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Antigen affinity purification.
Applications	IHC
Species Reactivity	Hu
Specificity	The antibody detects endogenous levels of total FOXR2 protein.
Immunogen Type	Recombinant Protein
Immunogen Description	Full length fusion protein
Target Name	FOXR2
Other Names	FOXN6
Accession No.	Swiss-Prot#: Q6PJJQ5NCBI Gene ID: 139628Gene Accssion: BC012934
Uniprot	Q6PJJQ5
GeneID	139628;
Concentration	0.9mg/ml
Formulation	Rabbit IgG in pH7.4 PBS, 0.05% NaN ₃ , 40% Glycerol.
Storage	Store at -20°C


Application Details

Immunohistochemistry: 1:25-1:100

Images



Immunohistochemical analysis of paraffin-embedded Human lung cancer tissue using #35740 at dilution 1/25.



Immunohistochemical analysis of paraffin-embedded Human thyroid cancer tissue using #35740 at dilution 1/25.

Background

Forkhead-box (FOX) transcription factors are implicated in carcinogenesis through gene amplification, retroviral integration, or chromosomal translocation. It expressed in breast cancer cell lines and primary cancer. The Subcellular location is in nucleus.

Note: This product is for in vitro research use only