

GST tag Polyclonal Antibody Cy5 Conjugated

Catalog No: #C03959Cy5

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	GST tag Polyclonal Antibody Cy5 Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Purified by Protein A.
Applications	IF
Species Reactivity	Ms GST fusion protein
Crossing Reactivity	GST tag fusion protein
Immunogen Description	Recombinant full length GST tag protein
Conjugates	Cy5
Target Name	GST tag
Other Names	glutathione-S-transferase; Glutathione S Transferase; Glutathione S transferase Mu 1; GST 26; GST class mu 1; GST1; GSTM1; GSTM1 1; GSTM1a 1a; GSTM1b 1b; GTH4; GTM1; H B; HB subunit 4; MGC26563; MU 1; MU; SJ26 antigen; SjGST.
Accession No.	NCBI Gene ID:2944
Uniprot	P09488
GeneID	2944;
Excitation Emission	625,650nm 670nm
Concentration	1mg ml
Formulation	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

Application Details

IF=1:50-200

Background

GST is a 26kDa protein encoded by the parasitic helminth *Schistosoma japonicum* and widely used in the pGEX family of GST plasmid expression vectors as a fusion protein with foreign proteins.

Note: This product is for in vitro research use only