

VBP1 antibody

Catalog No: #22061

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description


Product Name	VBP1 antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Purified by antigen-affinity chromatography.
Applications	WB
Species Reactivity	Hu
Immunogen Type	Recombinant protein
Immunogen Description	Recombinant protein fragment contain a sequence corresponding to a region within amino acids 1 and 143 of VBP1
Target Name	VBP1
Accession No.	Swiss-Prot:P61758Gene ID:7411
Uniprot	P61758
GeneID	7411;
Concentration	1mg/ml
Formulation	Supplied in 0.1M Tris-buffered saline with 20% Glycerol (pH7.0). 0.01% Thimerosal was added as a preservative.
Storage	Store at -20°C for long term preservation (recommended). Store at 4°C for short term use.

Application Details

Predicted MW: 23kd

Western blotting: 1:500-1:3000

Images



Sample (30 ug of whole cell lysate)
 A: A549
 12% SDS PAGE
 Primary antibody diluted at 1: 1000

Background

The protein encoded by this gene interacts with the Von Hippel-Lindau protein to form an intracellular complex. Because it functions as a chaperone protein, it is suspected that it may play a role in the transport of the Von Hippel-Lindau protein from the perinuclear granules to the nucleus or

cytoplasm. [provided by RefSeq]

Note: This product is for in vitro research use only