

DNAJB2 antibody

Catalog No: #22071

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description


Product Name	DNAJB2 antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Purified by antigen-affinity chromatography.
Applications	WB
Species Reactivity	Hu
Immunogen Type	Recombinant protein
Immunogen Description	Recombinant protein fragment contain a sequence corresponding to a region within amino acids 1 and 208 of DNAJB2
Target Name	DNAJB2
Other Names	HSJ1; HSPF3
Accession No.	Swiss-Prot:P25686Gene ID:3300
Uniprot	P25686
GeneID	3300;
Concentration	1mg/ml
Formulation	Supplied in 0.1M Tris-buffered saline with 10% Glycerol (pH7.0). 0.01% Thimerosal was added as a preservative.
Storage	Store at -20°C for long term preservation (recommended). Store at 4°C for short term use.

Application Details

Predicted MW: 36kd

Western blotting: 1:500-1:3000

Images



Sample (30 ug of whole cell lysate)

A: A431

10% SDS PAGE

DNAJB2 antibody diluted at 1: 1000

Note: This product is for in vitro research use only