

CD137 Polyclonal Antibody Cy3 Conjugated

Catalog No: #C04116Cy3

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	CD137 Polyclonal Antibody Cy3 Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Purified by Protein A.
Applications	IF
Species Reactivity	Hu Ms
Immunogen Description	KLH conjugated synthetic peptide aa 115-165 255 derived from human CD137
Conjugates	Cy3
Target Name	CD137
Other Names	ILA; 4-1BB; CD137; CDw137; Tumor necrosis factor receptor superfamily member 9; 4-1BB ligand receptor; T-cell antigen 4-1BB homolog; T-cell antigen ILA; TNFRSF9
Accession No.	Swiss-Prot#:Q07011NCBI Gene ID:3604
Uniprot	Q07011
GeneID	3604;
Excitation Emission	512,550nm 570,615nm
Cell Localization	Membrane bound, Cell membrane
Concentration	1mg ml
Formulation	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

Application Details

IF=1:50-200

Background

Receptor for TNFSF9 4-1BBL. Possibly active during T cell activation.

Note: This product is for in vitro research use only