

TEKT4 Antibody

Catalog No: #35960

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	TEKT4 Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Antigen affinity purification.
Applications	IHC
Species Reactivity	Hu
Specificity	The antibody detects endogenous levels of total TEKT4 protein.
Immunogen Type	Recombinant Protein
Immunogen Description	Fusion protein corresponding to residues near the C terminal of human tektin 4
Target Name	TEKT4
Other Names	tektin 4; tektin-4
Accession No.	Swiss-Prot#: Q8WW24NCBI Gene ID: 150483Gene Accssion: BC021716
Uniprot	Q8WW24
GeneID	150483;
Concentration	0.5mg/ml
Formulation	Rabbit IgG in pH7.4 PBS, 0.05% NaN ₃ , 40% Glycerol.
Storage	Store at -20°C

Application Details

Immunohistochemistry: 1:5-1:20

Images



Immunohistochemical analysis of paraffin-embedded Human lung cancer tissue using #35960 at dilution 1/10.

Background

Structural component of ciliary and flagellar microtubules. Forms filamentous polymers in the walls of ciliary and flagellar microtubules. Found in the abaxial (convex) surface of outer dense fibers in sperm flagella. TEKT4 is necessary for the proper coordinated beating of the sperm flagellum and male reproductive physiology.

Note: This product is for in vitro research use only