

C1QTNF6 Antibody

Catalog No: #35977

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	C1QTNF6 Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Antigen affinity purification.
Applications	IHC
Species Reactivity	Hu
Specificity	The antibody detects endogenous levels of total C1QTNF6 protein.
Immunogen Type	Recombinant Protein
Immunogen Description	Fusion protein corresponding to a region derived from internal residues of human C1q and tumor necrosis factor related protein 6
Target Name	C1QTNF6
Other Names	CTFP6; CTRP6; ZACRP6
Accession No.	Swiss-Prot#: Q9BXI9NCBI Gene ID: 114904Gene Accssion: BC020551
Uniprot	Q9BXI9
GeneID	114904;
Concentration	0.9mg/ml
Formulation	Rabbit IgG in pH7.4 PBS, 0.05% NaN ₃ , 40% Glycerol.
Storage	Store at -20°C

Application Details

Immunohistochemistry: 1:25-1:100

Images



Immunohistochemical analysis of paraffin-embedded Human brain tissue using #35977 at dilution 1/30.

Background

This protein is not uniformly distributed but rather focally localized in hepatocellular carcinoma cells and the tumor endothelial cells. It may play a role in metabolism and immunity. CTRPs acts as new molecular mediators connecting inflammatory and metabolic diseases.

Note: This product is for in vitro research use only