

IL-18R alpha CD218a Polyclonal Antibody Cy5 Conjugated

Catalog No: #C04222Cy5

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	IL-18R alpha CD218a Polyclonal Antibody Cy5 Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Purified by Protein A.
Applications	Flow-Cyt IF
Species Reactivity	HuB MsB RtB B
Immunogen Description	KLH conjugated synthetic peptide aa 23-73 541 derived from human IL-18R alpha
Conjugates	Cy5
Target Name	IL-18R alpha CD218a
Other Names	CD218a; IL18RA; IL1RRP; CDw218a; IL-1Rrp; Interleukin-18 receptor 1; IL-18R-1; IL-18R1; CD218 antigen-like family member A; IL1 receptor-related protein; IL1R-rp; IL18R1
Accession No.	Swiss-Prot#:Q13478NCBI Gene ID:8809
Uniprot	Q13478
GeneID	8809;
Excitation Emission	625,650nm 670nm
Cell Localization	Extracellular
Concentration	1mg ml
Formulation	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

Application Details

Flow-Cyt=1:50-200B IF=1:50-200B

Background

Receptor for interleukin 18 (IL-18). Binding to the agonist leads to the activation of NF-kappa-B.

Note: This product is for in vitro research use only