

# GBP5 antibody

Catalog No: #22085



Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)  
Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)


## Description

Product Name	GBP5 antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Purified by antigen-affinity chromatography.
Applications	WB
Species Reactivity	Hu
Immunogen Type	Recombinant protein
Immunogen Description	Recombinant protein fragment contain a sequence corresponding to a region within amino acids 1 and 183 of Human GBP5
Target Name	GBP5
Accession No.	Swiss-Prot:Q96PP8Gene ID:115362
Uniprot	Q96PP8
GeneID	115362;
Concentration	1mg/ml
Formulation	Supplied in 0.1M Tris-buffered saline with 10% Glycerol (pH7.0). 0.01% Thimerosal was added as a preservative.
Storage	Store at -20°C for long term preservation (recommended). Store at 4°C for short term use.

## Application Details

Predicted MW: 67kd  
Western blotting: 1:500-1:3000

## Images



Sample (30 ug of whole cell lysate)  
A: 293T  
7.5% SDS PAGE  
Primary antibody diluted at 1: 3000

Note: This product is for in vitro research use only