

NDUFS8 antibody

Catalog No: #22088



Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com


Description

Product Name	NDUFS8 antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Purified by antigen-affinity chromatography.
Applications	WB IHC
Species Reactivity	Hu
Immunogen Type	Recombinant protein
Immunogen Description	Recombinant protein fragment contain a sequence corresponding to a region within amino acids 1 and 208 of NDUFS8
Target Name	NDUFS8
Accession No.	Swiss-Prot:O00217Gene ID:4728
Uniprot	O00217
GeneID	4728;
Concentration	0.4mg/ml
Formulation	Supplied in 0.1M Tris-buffered saline with 20% Glycerol (pH7.0). 0.01% Thimerosal was added as a preservative.
Storage	Store at -20°C for long term preservation (recommended). Store at 4°C for short term use.


Application Details

Predicted MW: 24kd
Western blotting: 1:500-1:3000
Immunohistochemistry: 1:100-1:250

Images



Sample (30 ug of whole cell lysate)
A: Hep G2
12% SDS PAGE
NDUFS8 antibody diluted at 1: 1000



Immunohistochemical analysis of paraffin-embedded CA922 xenograft, using NDUFS8 antibody at 1: 250 dilution.

Note: This product is for in vitro research use only