

MC5-R antibody

Catalog No: #22126

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description


Product Name	MC5-R antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Purified by antigen-affinity chromatography.
Applications	WB
Species Reactivity	Hu
Immunogen Type	Peptide
Immunogen Description	Synthetic peptide contain a sequence corresponding to a region within amino acids 1 and 15 of MC5-R
Target Name	MC5-R
Accession No.	Swiss-Prot:P33032Gene ID:4161
Uniprot	P33032
GeneID	4161;
Concentration	0.2mg/ml
Formulation	Supplied in 0.1M Tris-buffered saline with 10% Glycerol (pH7.0). 0.01% Thimerosal was added as a preservative.
Storage	Store at -20°C for long term preservation (recommended). Store at 4°C for short term use.

Application Details

Predicted MW: 37kd

Western blotting: 1:500-1:3000

Images



Sample (30 ug of whole cell lysate)

A: HeLa

10% SDS PAGE

Primary antibody diluted at 1: 1000

Note: This product is for in vitro research use only