

CD64 IGFR1 Polyclonal Antibody Cy7 Conjugated

Catalog No: #C04753Cy7

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	CD64 IGFR1 Polyclonal Antibody Cy7 Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Purified by Protein A.
Applications	FCM IF(IHC-P)
Species Reactivity	Hu Ms Rt
Immunogen Description	KLH conjugated synthetic peptide aa 100-150 374 derived from human CD64 IGFR1
Conjugates	Cy7
Target Name	CD64 IGFR1
Other Names	CD64; FCRI; CD64A; IGFR1; High affinity immunoglobulin gamma Fc receptor I; IgG Fc receptor I; Fc-gamma RI; Fc-gamma RIA; FcgammaRIa; FCGR1A; FCG1; FCGR1
Accession No.	Swiss-Prot#:P12314NCBI Gene ID:2209
Uniprot	P12314
GeneID	2209;
Excitation Emission	743nm 767nm
Cell Localization	Extracellular
Concentration	1mg ml
Formulation	Aqueous buffered solution containing 1% BSA, 50% glycerol and 0.09% sodium azide.
Storage	Store at 4C for 12 months.

Application Details

FCM:1:20-100IF:1:50-200

Background

High affinity receptor for the Fc region of immunoglobulins gamma. Functions in both innate and adaptive immune responses.

Note: This product is for in vitro research use only