

PHEX antibody

Catalog No: #22133

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description


Product Name	PHEX antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Purified by antigen-affinity chromatography.
Applications	WB
Species Reactivity	Hu
Immunogen Type	Recombinant protein
Immunogen Description	Recombinant protein fragment contain a sequence corresponding to a region within amino acids 1 and 325 of Human PHEX
Target Name	PHEX
Accession No.	Swiss-Prot:P78562Gene ID:5251
Uniprot	P78562
GeneID	5251;
Concentration	1mg/ml
Formulation	Supplied in 0.1M Tris-buffered saline with 10% Glycerol (pH7.0). 0.01% Thimerosal was added as a preservative.
Storage	Store at -20°C for long term preservation (recommended). Store at 4°C for short term use.

Application Details

Predicted MW: 86kd

Western blotting: 1:500-1:3000

Images



Sample(30 ug whole cell lysate)
 A: MOLT4
 7.5% SDS PAGE
 Primary antibody diluted at 1: 1000

Background

The protein encoded by this gene is a transmembrane endopeptidase that belongs to the type II integral membrane zinc-dependent endopeptidase family. The protein is thought to be involved in bone and dentin mineralization and renal phosphate reabsorption. Mutations in this gene cause X-linked

hypophosphatemic rickets. [provided by RefSeq]

Note: This product is for in vitro research use only