

NGFRAP1 Antibody

Catalog No: #36629

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	NGFRAP1 Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Antigen affinity purification.
Applications	IHC
Species Reactivity	Hu
Specificity	The antibody detects endogenous levels of total NGFRAP1 protein.
Immunogen Type	Recombinant Protein
Immunogen Description	Full length fusion protein
Target Name	NGFRAP1
Other Names	Bex; BEX3; NADE; HGR74; DXS6984E
Accession No.	Swiss-Prot#: Q00994NCBI Gene ID: 27018Gene Accssion: BC003190/Q00994
Uniprot	Q00994
GeneID	27018;
Concentration	2.4mg/ml
Formulation	Rabbit IgG in pH7.4 PBS, 0.05% NaN ₃ , 40% Glycerol.
Storage	Store at -20°C


Application Details

Immunohistochemistry: 1:50-1:200

Images



Immunohistochemical analysis of paraffin-embedded Human brain tissue using #36629 at dilution 1/40.



Immunohistochemical analysis of paraffin-embedded Human thyroid cancer tissue using #36629 at dilution 1/40.

Background

Protein BEX3 is a protein that in humans is encoded by the NGFRAP1 gene. May be a signaling adapter molecule involved in p53-mediated apoptosis induced by NGF. Plays a role in zinc-triggered neuronal death. By similarity. May play an important role in the pathogenesis of neurodegenerative diseases.

Note: This product is for in vitro research use only