

## PITPNM1 Antibody

Catalog No: #36651

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)


## Description

|                       |                                                                                                                                    |
|-----------------------|------------------------------------------------------------------------------------------------------------------------------------|
| Product Name          | PITPNM1 Antibody                                                                                                                   |
| Host Species          | Rabbit                                                                                                                             |
| Clonality             | Polyclonal                                                                                                                         |
| Purification          | Antigen affinity purification.                                                                                                     |
| Applications          | IHC                                                                                                                                |
| Species Reactivity    | Hu                                                                                                                                 |
| Specificity           | The antibody detects endogenous levels of total PITPNM1 protein.                                                                   |
| Immunogen Type        | Recombinant Protein                                                                                                                |
| Immunogen Description | Fusion protein corresponding to residues near the C terminal of human phosphatidylinositol transfer protein, membrane-associated 1 |
| Target Name           | PITPNM1                                                                                                                            |
| Other Names           | Rd9; NIR2; RDGB; DRES9; RDGB1; RDGBA; PITPNM; RDGBA1                                                                               |
| Accession No.         | Swiss-Prot#: O00562NCBI Gene ID: 9600Gene Accssion: BC022230/O00562                                                                |
| Uniprot               | O00562                                                                                                                             |
| GeneID                | 9600;                                                                                                                              |
| Concentration         | 2mg/ml                                                                                                                             |
| Formulation           | Rabbit IgG in pH7.4 PBS, 0.05% NaN <sub>3</sub> , 40% Glycerol.                                                                    |
| Storage               | Store at -20°C                                                                                                                     |


## Application Details

Immunohistochemistry: 1:25-1:100

## Images



Immunohistochemical analysis of paraffin-embedded Human thyroid cancer tissue using #36651 at dilution 1/40.



Immunohistochemical analysis of paraffin-embedded Human cervical cancer tissue using #36651 at dilution 1/40.

## Background

PITPNM1 belongs to a family of membrane-associated phosphatidylinositol transfer domain-containing proteins that share homology with the *Drosophila* retinal degeneration B (rdgB) protein.

Note: This product is for in vitro research use only