

LOC169436 antibody

Catalog No: #22146



Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com


Description

Product Name	LOC169436 antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Purified by antigen-affinity chromatography.
Applications	WB IF
Species Reactivity	Hu
Immunogen Type	Recombinant protein
Immunogen Description	Recombinant protein fragment contain a sequence corresponding to a region within amino acids 438 and 624 of Human C9orf96
Target Name	LOC169436
Accession No.	Swiss-Prot:Q8NE28Gene ID:169436
Uniprot	Q8NE28
GeneID	169436;
Concentration	1mg/ml
Formulation	Supplied in 0.1M Tris-buffered saline with 10% Glycerol (pH7.0). 0.01% Thimerosal was added as a preservative.
Storage	Store at -20°C for long term preservation (recommended). Store at 4°C for short term use.


Application Details

Predicted MW: 76kd
Western blotting: 1:500-1:3000
Immunofluorescence: 1:100-1:200


Images



Sample (30 ug of whole cell lysate)
A: A431
7.5% SDS PAGE
Primary antibody diluted at 1: 5000



Merged with DNA probe



Immunofluorescence analysis of paraformaldehyde-fixed HeLa, using LOC169436 antibody at 1: 200 dilution.

Note: This product is for in vitro research use only