



Introduction 0

Technical Support 0

1 General Biological Materials

PCR Enzymes 1

- Taq DNA Polymerases
- RT Enzymes
- qPCR Products

SafeView™ Products 4

DNA Marker & Protein Markers 4

Western Blot and ELISA Products 5

Transfection Reagents 5

Other Enzymes Used in Molecular Biology..... 6

Protein Purification Products 7

2 Cellular Materials

Immortalized Cells 9

Primary Cells 10

Cell Culture Related Products 10

Mycoplasma Detection Kit 11

Cell Immobilization Agents 11

Stem Cells 12

Antibodies 12

- Gene-Specific Antibodies
- Tag Antibodies
- Loading Control Antibodies

Recombinant Proteins and Protein Lysates (Human, Mouse, Rat) 14

3 Genetic Materials

ORF cDNA Library (Human, Mouse, Rat).....	15
ORF Lentiviral Vector Library(Human, Mouse, Rat)	15
ORF Adenoviral Library (Human, Mouse, Rat)	15
ORF Cumate inducible Lentiviral Library (Human, Mouse, Rat)	15
iLenti siRNA Library (Human, Mouse, Rat)	16
miRNA Library (Human, Mouse, Rat)	16
ORF Protein Expression Vector Library(Human, Mouse, Rat)	16

4 Custom Services

Custom Cloning	17
Site-direct Mutagenesis	17
Lentivirus	17
Adenovirus	17
Stable Cell Line	17
Custom Antibody	17
Recombinant Protein production	17
Custom Cell Immortalization	18
miRNA Profiling	18
Whole Genome Sequencing	18
Plasmid Sequencing	18
Oligo Synthesis	18

Introduction

abm is a world leader in developing and marketing research reagents for both academic and biotech communities. Our company is evolving into an one-stop solution for any biological material you may encounter in any project. We are currently the only whole-sale discount supplier in biological materials and are happy to be able to provide you with top value general, cellular, and genetic materials. As our customer, you can rely on us to save you over 50% on any biological material without any compromise on product quality. In addition to the discounts we offer, we are also very innovative in providing thousands of unique biological materials that no other companies currently provide, and with our continual and rapid growth many more novel products will be added to our product line in the years ahead. When you require any material for your research project, whether a gene or protein, simply explore our comprehensive website (www.abmGood.com) to find the best value and innovative products to expedite your ongoing project.

Another great advantage of becoming a valued **abm** customer is that every product **abm** offers will be shipped with a rewarding gift as appreciation of your customer. Please contact your local sale representatives or our customer service associates for more details.

Technical Support

Applied Biological Materials Inc.

Email Addresses:



General Information: info@abmGood.com

Order Products: order@abmGood.com

Technical Support: technical@abmGood.com

Business Development: bd@abmGood.com

Telephone: (8:30am - 4:30pm PST, Mon - Fri)



Toll Free: 1-866-757-2414

Local: 604-247-2416

Fax: 604-247-2414 (24 hours)

Address:

Suite #8-13520 Crestwood Place,
Richmond, BC, Canada V6V 2G2

Website: www.abmGood.com

1. General Biological Materials

abm's general biological products cover all materials used daily in any lab for genetic and cellular projects. This includes enzymes for PCR, RT-PCR, qPCR, and other recombinant DNA technologies, innovative DNA staining dyes (SafeView Plus™), markers for DNA and protein, transfection reagents, and small chemical molecules.

1.1 PCR Enzymes

abm is the leader in next generation PCR enzymes for both quality and discount. All of our enzymes are designed with next generation enzyme technology for highest sensitivity, fidelity, and specificity. They are the premier choices for whole genome sequencing and PCR assays used in personalized medicine.

1.1.1 Taq DNA Polymerases

The most comprehensive range of polymerases for any PCR application, including wild-type Taq Polymerase, and different mutant forms of Taq polymerase for fast PCR, high fidelity, long templates and dirty templates. The availability of our PCR-Sure™ Kits will guarantee universal PCR success for even the most challenging templates.

Product Name	Size	Cat. No.	Price
Bestaq DNA Polymerase	100 U (50 µl)	G456	\$65.00
Taq DNA Polymerase	1000 U (200 µl)	G009	\$40.00
Taq Plus DNA Polymerase	250 U (50 µl)	G012	\$35.00
Precision™ Taq DNA Polymerase	500 U (100 µl)	G078	\$150.00
HotStart Taq DNA Polymerase	250 U (50 µl)	G011	\$62.50
2X PCR Bestaq MasterMix	5 ml (200 PCR reactions)	G464	\$185.00
2X PCR Taq MasterMix	5 ml (200 PCR reactions)	G013	\$115.00
2X PCR Taq Plus MasterMix	5 ml (200 PCR reactions)	G014	\$135.00
TaqFast DNA Polymerase PCR	250 U (50 µl)	G277	\$62.50
PCR-Sure™ Kit	12 x 10 reactions	G065	\$169.00
Bloodirect DNA Polymerase	100 U (100 µl)	G462	\$40.00

1.1.2 RT Enzymes

abm provides Advanced, improved functional Reverse Transcriptases (RT) for every application to meet your cDNA synthesis needs, including long and complicated RNA templates with secondary structure.

Product Name	Size	Cat. No.	Price
EasyScript™ RT	25 reactions	G231	\$35.00
EasyScript Plus™ RT	25 reactions	G177	\$65.00
EasyScript™ cDNA Synthesis Kit	25 reactions	G233	\$65.00
EasyScript Plus™ cDNA Synthesis Kit	25 reactions	G235	\$90.00
EasyScript™ cDNA Synthesis SuperMix	25 reactions	G451	\$60.00
EasyScript Plus™ cDNA Synthesis SuperMix	25 reactions	G453	\$85.00

1.1.3 qPCR Products

Sensitive, reliable qPCR MasterMixes for any brand of qPCR equipment at discount prices. All qPCR products are offered as 500 reactions size (4 x 1.25ml). EvaGreen MasterMixes are used for qPCR reactions based on DNA-binding fluorescent dye (ex. SYBR Green I); TaqProbe MasterMixes are used for qPCR reaction based on TaqMan probe; Whereas KiloGreen MasterMix are used for qPCR reactions with long amplicon up to 1.3kb, therefore could be used in the library construction for whole genome sequencing (WGS).

Cat. No.	Product Name	qPCR Instruments	Price
MasterMix-R	EvaGreen 2X qPCR MasterMix-ROX	-ABI® 7000, 7300, 7700, 7900, 7900HT, StepOnePlus™, StepOne™, OpenArray, PRISM™ Sequencing Detection Series	\$155.00
MasterMix-P	TaqProbe 2X qPCR MasterMix-ROX	-Biometra Topical -Fluidigm BioMark™	\$155.00
MasterMix-KR	KiloGreen 2X qPCR MasterMix-ROX	-Watogene SmartChip System -TianLong TL998 System	\$195.00
MasterMix-LR	EvaGreen 2X qPCR MasterMix-Low ROX	-ABI® 7500 ViiA™, QuantStudio	\$155.00
MasterMix- PL	TaqProbe 2X qPCR MasterMix -Low ROX	-BioGene InSyte™ -Illumina Eco	\$155.00
MasterMix-KL	KiloGreen 2X qPCR MasterMix-Low ROX	-Stratagene® Mx3000, Mx3005, Mx4000 -Analytikjena qTower Series	\$195.00
MasterMix-iC	EvaGreen 2X qPCR MasterMix-iCycler		\$155.00
MasterMix-PC	TaqProbe 2X qPCR MasterMix-iCycler	-BioRad® iCycler®, iQ™5, MyiQ™	\$155.00
MasterMix-KC	KiloGreen 2X qPCR MasterMix-iCycler		\$195.00
MasterMix-S	EvaGreen 2X qPCR MasterMix	-BioRad® CFX96, CFX384, Chromo4™, CFX Connect™, Opticon 2, MiniOpticon™ -Roche LightCycler® (2.0, 1.5, 480, 1536, Nano) -MJ Research Opticon™, Opticon™ 2, Chromo® 4	\$155.00
MasterMix-PS	TaqProbe 2X qPCR MasterMix-no dye	-Eppendorf® Realplex 4 -BioGene SynChron™ -Corbett Rotor-gene® (3000, 6200, 62H0, 6500, 65H0, 6600) -Eppendorf Mastercycler® realplex (s, 4, 4s), Pro (S, 384), Nexus (gradient, eco, flat) -Cepheid SmartCycler®, GeneXpert	\$155.00
MasterMix-KS	KiloGreen 2X qPCR MasterMix	-Enigma® ML -Idaho LightScanner® (24, 32), RapidCycler®2, R.A.P.I.D (LT, LT Food), RAZOR EX, JBAIDS -Qiagen Rotor-Gene™ (Q, 6000) -Takara Dice™ -Thermo Scientific PikoReal -DNA-Technology DT96, DTlite, DT-322 -Bioer LineGene (3310/3320, K FQD-48A, I, II, 9620, 9640, 9660, 9680) -Bioneer Exicycler™	\$195.00
MasterMix-PM	TaqProbe 2X qPCR MasterMix-Multiplex	Any qPCR instrument that supports multiplex reactions	\$155.00

1.2 SafeView™ products

Carcinogenic toxicity associated with Ethidium Bromide gel electrophoresis is a major health and environmental concern. Scientists at **abm** have successfully developed alternative, sensitive, non carcinogenic SafeView™ DNA/RNA staining dyes of different colors (SafeView™, SafeView Plus™, Safe-Green™, Safe-Red™, and Safe-White™). SafeView Plus™ is the most advanced and safe DNA staining dye, it has been proven to be superior to any other similar DNA staining dyes on the market.

Product Name	Size	Cat. No.	Price
SafeView Plus™	1 ml (10,000X)	G468	\$60.00
Safe-Pack™	G108 - W,G,R	G108-P	\$55.00
SafeView™	1 ml	G108	\$25.00
Safe-Green™	1 ml	G108-G	\$25.00
Safe-Red™	1 ml	G108-R	\$25.00
Safe-White™	1 ml	G108-W	\$25.00

1.3 DNA Marker & Protein Markers

Product Name	Range	Bands	Size	Cat. No.	Price
100 bp Opti-DNA Marker	50 bp - 1.5 kb	17	500 µl/100 loads	G016	\$75.00
100 bp Plus Opti-DNA Marker	100 bp - 3 kb	12	500 µl/100 loads	G193	\$75.00
1 kb Opti-DNA Marker	100 bp - 10 kb	19	500 µl/100 loads	G106	\$75.00
1 kb Plus Opti-DNA Marker	250 bp - 25 kb	14	500 µl/100 loads	G248	\$75.00
Opti-Protein Marker	10 kDa - 175 kDa	11	500 µl	G252	\$95.00
Opti-Protein XL Marker	10 kDa - 245 kDa	12	500 µl	G266	\$110.00

1.4 Western Blot and ELISA Products

Quality, discount, whole-sale supplier for ECL, ELISA substrates, and hundreds of secondary antibodies for both Western blot and ELISA applications.

Product Name	Size	Cat. No.	Price
1X Protein Loading Buffer	5 ml	G032	\$35.00
5X Protein Loading Buffer	3 ml	G031	\$45.00
ECL Western Blotting Detection Reagents	2 X 100 ml	G075	\$90.00
Protease Inhibition Cocktail with EDTA	1 ml	G135E	\$100.00
Protease Inhibitor Cocktail	1 ml	G135	\$70.00
ELISA ABTS Substrate for HRP	250 ml	G445	\$75.00
ELISA TMB Substrate for HRP	250 ml	G443	\$95.00
ELISA Ultrasensitive TMB Substrate for HRP	250 ml	G444	\$85.00
Horseradish Peroxidase	100 mg	G015	\$100.00
BCIP-NBT Dye for IHC	250 ml	G447	\$175.00
Fast Red Violet Dye for IHC	250 ml	G446	\$250.00
DTT	3 X 1 ml	G025	\$35.00
PMSF	100 mM in isopropanol; 3 X 1 ml	G034	\$35.00

1.5 Transfection Reagents

Premium transfection reagents offer the highest transduction efficiency and minimal toxicity.

Product Name	Size	Cat. No.	Price
DNAfectin™ 2100 Transfection Reagent	1 mg/ml	G2100	\$145.00
RNAfectin™ Transfection Reagent	1 ml (100 - 500 transfections)	G073	\$165.00
Lentifectin™ Transfection Reagent	1 ml	G074	\$155.00
Calciumfectin™ Mammalian Transfection Reagent Kit	Kit (100 transfections)	G099	\$125.00
Profectin™ Transfection Reagent	250 µl	G288	\$220.00

1.6 Other Enzymes Used in Molecular Biology

Quality enzymes have been used in whole genome sequencing and general recombinant DNA technologies for performance and significantly discount prices.

Next Generation Sequencing (NGS) Enzymes and Tools

Product Name	Size	Cat. No.	Price
<i>E.coli</i> DNA ligase	100ul (1000U, 10U/ul)	E007	\$185
T4 DNA Ligase	100ul (100U, 1 U/ul)	G467	\$55
T4 RNA Ligase 1 (ssRNA Ligase)	100ul (1000U, 10U/ul)	E010	\$35
T4 RNA Ligase 2 (truncated)	100ul (20000U, 200U/ul)	E012	\$325
<i>E.coli</i> DNA Polymerase I	100ul (1000U, 10U/ul)	G469	\$95
Klenow Fragment (3'→5' exo-)	100ul (500U, 5U/ul)	E009	\$95
<i>E.coli</i> DNA Polymerase I, Large (Klenow) Fragment	100ul (500U, 5U/ul)	E013	\$95
Phi29 DNA Polymerase	100ul (1000U, 10U/ul)	E014	\$225
T4 DNA Polymerase	100ul (300U, 3U/ul)	E005	\$75
RNase H, <i>E.coli</i>	100ul (500U, 5U/ul)	E018	\$85
RNase Inhibitor	100ul (4000U, 40U/ul)	G138	\$75

Cloning kits

Product Name	Size	Cat. No.	Price
Ligation-Free Cloning System	25 reactions	E001	\$295
Ligation-Free Cloning System	100 reactions	E002	\$980
ProClone™ Competent Cell	5ml	E003	\$95
Cloning Optimizer	200ul	E004	\$95

DNA/RNA Modifying Enzymes

Product Name	Size	Cat. No.	Price
DNA/RNA Ligases			
<i>E.coli</i> DNA ligase	100ul (1000U, 10U/ul)	E007	\$185
T4 DNA Ligase	100ul (100U, 1 U/ul)	G467	\$55
T4 RNA Ligase 1 (ssRNA Ligase)	100ul (1000U, 10U/ul)	E010	\$35
T4 RNA Ligase 2 (dsRNA Ligase)	100ul (1000U, 10U/ul)	E011	\$225
T4 RNA Ligase 2 (truncated)	100ul (20000U, 200U/ul)	E012	\$325
Mesophilic DNA Polymerases			
<i>E.coli</i> DNA Polymerase I	100ul (1000U, 10U/ul)	G469	\$95
<i>E.coli</i> DNA Polymerase I, Large (Klenow) Fragment	100ul (500U, 5U/ul)	E013	\$95
Klenow Fragment (3'→5' exo-)	100ul (500U, 5U/ul)	E009	\$95
Phi29 DNA Polymerase	100ul (1000U, 10U/ul)	E014	\$225
T4 DNA Polymerase	100ul (300U, 3U/ul)	E005	\$75
Nucleases			
Exonuclease I, <i>E.coli</i>	100ul (2000U, 20U/ul)	E015	\$45
Exonuclease III, <i>E.coli</i>	100ul (10000U, 100U/ul)	E016	\$85
RNases			
RNase H, <i>E.coli</i>	100ul (500U, 5U/ul)	E018	\$85
RNase I, <i>E.coli</i>	100ul (5000U, 50U/ul)	E019	\$45
RNase III, <i>E.coli</i>	100ul (200U, 2U/ul)	E020	\$85
Thermostable RNase H	100ul (500U, 5U/ul)	E008	\$145
RNase A	25mg	G117	\$20
Phosphatases and Sulfurylases			
Inorganic Pyrophosphatase, Thermostable	100ul (200U, 2U/ul)	E021	\$55
Kinases			
T4 Polynucleotide Kinase	100ul (1000U, 10U/ul)	E006	\$85
Single-Stranded DNA Binding Proteins			
Extreme Thermostable SSB	100ul, 1mg/ml	E023	\$125
<i>E.coli</i> SSB	100ul, 1mg/ml	E022	\$95
RecA Protein, <i>E.coli</i>	100ul, 1mg/ml	E024	\$35
RecA Protein, Tth	100ul, 1mg/ml	E025	\$95
Others			
RNase Inhibitor	100ul (4000U, 40U/ul)	G138	\$75
DpnI	100ul (2000U, 20U/ul)	E026	\$85
DNase I	1ml (2000U, 2U/ul)	G028	\$55

1.7 Protein Purification Products

We offer a complete range of different resins, His and Flag tag antibodies at most competitive prices.

Product Name	Size	Cat. No.	Price
Ni-IDA Agarose Resins for Protein Purification	5 ml	G250	\$55.00
Glutathione Agarose Bead Purification Resins	10 ml	G241	\$165.00
Protein A Agarose Column for antibody purification	1 X 1 ml	G245	\$125.00
Enterokinase	100 Units	G699	\$35.00

2. Cellular Materials

abm is globally recognized leader in primary and immortalized cell lines. We have the largest collection of immortalized human, mouse, and rat cell lines in the world. This includes immortalized human, mouse, rat hepatocytes that can't be found anywhere else. In addition, we also offer reliable and affordable services for any primary cell immortalization project and iPSC project. Our other cutting-edge cell culture products include stem cells (iPSC), viral vectors, mycoplasma detection kits, and general cell culture items such as serum and medium.

Within our cellular materials product range are a unique collection of Tag antibodies and Western Blot loading control antibodies. Because of their highest quality and discount, those products are now used by 70% worldwide customers. Besides these popular antibodies, we also have over 30,000 gene-specific antibodies and the list is growing every day. Another unique collection of cellular materials is recombinant proteins for any human, mouse, and rat gene. Our recombinant proteins are so affordable that there is now no reason for you to invest your precious time and further expense on any unpredictable protein purification project.

2.1 Immortalized Cells

Product Name	Size	Cat. No.	Price
Immortalized Human Carotid Artery Smooth Muscle Cells-hTERT	1x10 ⁶ cells/ml	T0511	Inquire
Immortalized Human Chondrocytes-SV40	1X10 ⁶ cells/ml	T0021	Inquire
Immortalized Human Microglia-SV40	1X10 ⁶ cells/ml	T0251	Inquire
Immortalized Hair Follicle Dermal Papilla Cells - SV40	1X10 ⁶ cells/ml	T0500	Inquire
Immortalized Human Cardiomyocytes - SV40	1x10 ⁶ cells/ml	T0445	Inquire
Immortalized Human Hepatocytes-SV40	1X10 ⁶ cells/ml	T0050	Inquire
Immortalized Human Diploid Fibroblasts -hTERT	1X10 ⁶ cells/ml	T1013	Inquire
Immortalized Human Ovarian Epithelial Cells-SV40	1X10 ⁶ cells/ml	T1074	Inquire
Immortalized Human Pancreas-Beta Cells-SV40	1X10 ⁶ cells/ml	T0170A	Inquire
Immortalized Human Benign Prostate Epithelial Cell-SV40	1X10 ⁶ cells/ml	T0002	Inquire
Immortalized Human Skin Fibroblast-Myc	1X10 ⁶ cells/ml	T0301	Inquire
Immortalized Human Vascular Endothelial Cells [EC-RF24]	1X10 ⁶ cells/ml	T0003	Inquire

2.2 Primary Cells

Product Name	Size	Cat. No.	Price
Human Primary Subcutaneous Adipocytes	1 X 10 ⁶ cells/ml	T4000	\$625.00
Human Primary Alveolar Cells-Type 1	1 X 10 ⁶ cells/ml	T4003	\$625.00
Peripheral Blood Mononuclear Cells	75 X 10 ⁶ cells/ml	T4103	\$525.00
Human Primary Aortic Endothelial Cells	1 X 10 ⁶ cells/ml	T4039	\$625.00
Human Primary Osteoblasts	1 X 10 ⁶ cells/ml	T4015	\$625.00
Human Primary Dendritic Cells	1 X 10 ⁶ cells/ml	T4021	\$550.00
Human Primary Astrocytes	1 X 10 ⁶ cells/ml	T4033	\$625.00
Human Primary Mammary Epithelial Cells	1 X 10 ⁶ cells/ml	T4036	\$625.00
Human Primary Colon Cells	1 X 10 ⁶ cells/ml	T4056	\$625.00
Human Primary Retinal Pigment Epithelial Cells	1 X 10 ⁶ cells/ml	T4057	\$625.00
Human Primary Hair Follicle Dermal Papilla Cells	1 X 10 ⁶ cells/ml	T4059	\$625.00
Human Primary Cardiomyocytes	1 X 10 ⁶ cells/ml	T4037	\$625.00
Human Primary Cryopreserved Non-Parenchymal Cells	1 X 10 ⁶ cells/ml	T4061	\$625.00
Human Primary Mesothelial Cells	1 X 10 ⁶ cells/ml	T4067	\$625.00
Human Primary Oral Epithelial Cells	1 X 10 ⁶ cells/ml	T4074	\$625.00
Human Primary Pancreatic Cells	1 X 10 ⁶ cells/ml	T0199	\$625.00
Human Primary Nasopharyngeal Cells	1 X 10 ⁶ cells/ml	T4070	\$625.00
Human Primary Placental Epithelial Cells	1 X 10 ⁶ cells/ml	T4076	\$625.00
Human Primary Mesangial Cells	1 X 10 ⁶ cells/ml	T4086	\$625.00
Human Primary Prostate Cells	1 X 10 ⁶ cells/ml	T4077	\$625.00
Human Primary Skeletal Muscle Myoblasts	1 X 10 ⁶ cells/ml	T4068	\$625.00
Human Primary Dermal Endothelial cells	1 X 10 ⁶ cells/ml	T4097	\$625.00
Human Primary Umbilical Vein Endothelial Cells	1 X 10 ⁶ cells/ml	T4098	\$625.00

2.3 Cell Culture Related Products

Product Name	Size	Cat. No.	Price
T25 Flasks	10 Flasks	G299	\$50.00
Prigrow I Medium	500 ml	TM001	\$95.00
Prigrow II Medium	500 ml	TM002	\$95.00
Prigrow III Medium	500 ml	TM003	\$95.00
Prigrow IV Medium	500 ml	TM004	\$95.00
Trypsin-EDTA	100 ml	TM051	\$15.00
Trypsin-EDTA	500 ml	TM050	\$25.00
USDA Research Grade Origin Fetal Bovine Serum	100 ml	TM999-100	\$55.00
Penicillin/Streptomycin	100 ml	G255	\$25.00

2.4 Mycoplasma Detection Kit

Product Name	Size	Cat. No.	Price
PCR Mycoplasma Detection Kit	100 Reactions	G238	\$175.00
Mycoplasma Elimination Cocktail	2 x 1.0ml	G398	\$175.00

2.5 Cell Immobilization Agents

Product Name	Description	Size	Cat. No.	Price
Lenti-SV40 T Antigen Virus	Recombinant Lentivirus containing SV40 large T antigen (10^6 cfu/ml)	10 ml	G256	\$750.00
Adeno-SV40 Virus	Recombinant Adenovirus containing SV40 large and small T antigen (10^6 pfu/ml)	250 μ l	G210	\$950.00
Retro-SV40 Virus	Recombinant Retrovirus containing SV40 large and small T antigen (10^6 cfu/ml)	10 ml	G212	\$850.00
Lenti-hTERT Virus	Recombinant hTERT Lentivirus (10^6 cfu/ml)	10 ml	G200	\$750.00
Adeno-hTERT	Recombinant Telomerase Adenovirus (10^6 pfu/ml)	250 μ l	G205	\$950.00
Retro-E1/hTERT Virus	Recombinant hTERT Retrovirus (10^6 cfu/ml)	10 ml	G207	\$850.00
Lenti-Myc T58A Virus	Recombinant Myc T58A Lentivirus (10^6 cfu/ml)	10 ml	G217	\$750.00
Lenti-Ras V12 Virus	Recombinant Ras V12 Lentivirus (10^6 cfu/ml)	10 ml	G221	\$750.00
Lenti-p53 siRNA Virus	Recombinant p53 siRNA Lentivirus (10^6 cfu/ml)	10 ml	G219	\$750.00
Lenti-Rb siRNA Virus	Recombinant Rb siRNA Lentivirus (10^6 cfu/ml)	10 ml	G223	\$750.00

2.6 Stem Cells

Product Name	Description	Cat. No.	Size	Price
iPSC-Minicircle DNA (Human)	Minicircle DNA-Yamanaka	G388	100 µg	\$330.00
	Minicircle DNA-Thomson	G389		
iPSC Lentivirus (Human)	EF1α Promoter Thomson Factors Lentivirus Set	G356	5 X 200 µl	\$775.00
	EF1α Promoter Yamanaka Factors Lentivirus Set	G361	1 X 10 ⁸ cfu/ml	
iPSC Adenovirus (Human)	Thomson Factor Set Adenovirus	000777A	5 X 200 µl	\$1,250.00
	Yamanaka Factor Set Adenovirus	000779A		
iPSC-Minicircle DNA (Mouse)	Minicircle DNA-mYamanaka	G390	100 µg	\$330.00
	Minicircle DNA-mYamanaka	G391		
iPSC Lentivirus (Mouse)	EF1α Promoter mThomson Factors Lentivirus Set	G371	5 X 200 µl	\$775.00
	EF1α Promoter mYamanaka Factors Lentivirus Set	G372	1 X 10 ⁸ cfu/ml	
iPSC Adenovirus (Mouse)	mThomson Factor Set Adenovirus	000787A	5 X 250 µl	\$3,250.00
	mYamanaka Factor Set Adenovirus	000789A		
CB iPSC Plasmids	pEV-SFFV-Oct4-2A-Sox2 Vector	G505	1 µg	\$275.00
	Lenti-SFFV-Oct4-2A-Sox2 Virus	G438	1 X 10 ⁶ cfu/ml	\$50.00

2.7 Antibodies

2.7.1 Gene-Specific Antibodies

Besides popular tag and loading control antibodies, **abm** also offers 30,000+ gene specific antibodies which is one of the largest in the world. Because of the largest ORF vector library we have, scientists at **abm** are in unique position to evaluate each antibody's specificity with lysates generated from 293 cells transfected with a specific ORF vector. Please visit our website www.abm-Good.com to search any specific antibody you are looking for. If you can't find the antibody you are interested in, please contact our customer service associates and you will be assisted instantly.

2.7.2 Tag Antibodies

Product Name	Size	Cat. No.	Price
His Tag Antibody(Mouse Monoclonal)	100 µg	G020	\$115.00
HA Tag Antibody(Mouse Monoclonal)	100 µg	G036	\$115.00
GFP Tag Antibody (Mouse Monoclonal)	100 µg	G096	\$115.00
GST Tag Antibody (Mouse Monoclonal)	100 µg	G018	\$115.00
Myc Tag Antibody (Mouse Monoclonal)	100 µg	G019	\$115.00
V5 Tag Antibody (Rabbit Polyclonal)	100 µl	G466	\$115.00
V5 Tag Antibody(Mouse Monoclonal)	100 µg	G189	\$115.00
D Tag Antibody (Mouse Monoclonal)	100 µg	G191	\$115.00
RFP Tag Antibody (Mouse Monoclonal)	100 µg	G093	\$195.00
MBP (Maltose Binding Protein) Tag Antibody	100 µg	G079	\$145.00
VSVG Tag Antibody (Goat Polyclonal)	100 µg	G192	\$255.00
YFP Tag Antibody (Mouse Monoclonal)	100 µg	G163	\$115.00
CFP Tag Antibody (Monoclonal)	100 µg	G164	\$115.00

2.7.3 Loading Control Antibodies

Product Name	Size	Cat. No.	Price
Anti-GAPDH Loading Control Monoclonal Antibody	100 µg	G041	\$115.00
Anti-GAPDH Polyclonal Loading Control Antibody	100 µl	Y058203	\$115.00
Anti-Beta-Actin Monoclonal Loading Control Antibody	100 µg	G043	\$115.00
Anti-Beta-Actin Polyclonal Loading Control Antibody	100 µg	G046	\$115.00
Anti-Alpha-Tubulin Monoclonal Loading Control Antibody	100 µg	G094	\$115.00
Anti-Beta-Tubulin Monoclonal Loading Control Antibody	100 µg	G098	\$115.00

2.8 Recombinant Proteins and Protein Lysates

(Human, Mouse, Rat)

Our extensive collection of ORF library also extend to protein vector. We have each ORF cDNA already cloned into proper protein expression vectors in a ready-to-express format. As a valuable customer, you have the choices of purchasing a protein expression vector with any ORF cDNA of your interest or expressed protein or cell lysates directly from us. Please visit our website www.abmGood.com to search any ORF protein expression vector, recombinant proteins, and positive lysates you are interested in. If you can't find your materials due to any reason, please contact our customer service associate and you will be assisted instantly.

3. Genetic Materials

abm provides the most comprehensive selection of genetic materials in the world, whether it is a cDNA from any specie, siRNA, miRNA, promoter, or in any form vectors of either viral or non-viral origin. We are the one stop solution for your needs of any genetic material. If any material you are looking for is not currently in our database, we can custom-generate it for you at material cost. We are the expert on recombinant DNA technologies and can do better than other competitors in the world.

3.1 ORF cDNA Library (Human, Mouse, Rat)

Open Reading Frame of cDNAs from any specie is within our product and service ranges. Please visit our website www.abmGood.com to search any ORF cDNA material you require. If you are unable to find your required materials for any reason, please contact our customer service associates.

3.2 ORF Lentiviral Vector Library (Human, Mouse, Rat)

We have a complete collection of lentivirus ORFs for any specie in our inventory, please visit our website www.abmGood.com to search any ORF cDNA material you are looking for. If you can't find your materials due to any reason, please contact our customer service associates and you will be assisted instantly.

3.3 ORF Adenoviral Library (Human, Mouse, Rat)

We have the world's only complete collection of adenovirus ORF for any specie in our inventory at **abm**. Please visit our website www.abmGood.com to search any ORF cDNA material you are looking for. If you can't find your materials for any reason, please contact our customer service associates and you will be assisted instantly.

3.4 ORF Cumate Inducible Lentiviral Library (Human, Mouse, Rat)

The cumate inducible expression system is the most advanced technology for regulated gene expression developed in the last few years by leading scientists at National Research Laboratories (NRC) of Canada. Due to our leading position in viral vectors, the NRC agreed to license us the right to this technology for world-wide commercialization in recombinant viral vectors. The system has been proven to be the most reliable with the tightest control

on background leak expression in comparison to other traditional inducible expression systems on the market. To date, human, mouse, and rat ORF clones of Cumate inducible Lentivirus have been generated for immediate applications as long as ORF inserts are less than 3.0 kb for a single vector format and 5.0 kb for a duo vector format. Inserts over the specified sizes will make the recombinant lentiviral vectors beyond the lentiviral vector packaging limit. Please visit our website www.abmGood.com for your specific inducible vector needs.

3.5 iLenti siRNA Library (Human, Mouse, Rat)

Scientists at **abm** have developed the most advanced, convergent promoter, signature iLenti siRNA for any human, mouse, rat genes. Simply visit our website www.abmGood.com to search for any siRNA lentiviral vector you may require. Advantages in using our convergent promoter system over a traditional hairpin loop single vector include:

- a. Flexibility for either transfection or viral infection.
- b. Vector stability no special competent cells required.
- c. Potent 27 - 29 bp oligos.
- d. Convergent promoters to avoid hairpin loop structures, making sequencing and plasmid growth much easier.

3.6 miRNA Library (Human, Mouse, Rat)

abm is one of few pioneering institutions that has developed a whole genome human and mouse miRNA lentivirus library and a qPCR-based whole genome miRNA profiling system. In addition, we have also developed the only library in the world for every human and mouse miRNA adenoviral vector. So, for any miRNA lentiviral and adenoviral vector as well as qPCR miRNA array, **abm** stands alone as the top choice for your project.

3.7 ORF Protein Expression Vector Library (Human, Mouse, Rat)

Because of our extensive collection for human, mouse, and rat ORF cDNAs, we have generated bacteria protein expression vectors for every ORF material we have. Please visit our website www.abmGood.com to search any ORF expression vector you require for immediate recombinant protein induction and purification.

4. Custom Services

Over the last 10 years, scientists at **abm** have developed a variety of special technology platforms and served thousands of customers for their special material needs in the following areas. The services we provided are so affordable (material and labour costs) and reliable that there is really no reason that you would like to take any of the following work yourselves. In addition, our superb services are risk-free and you will not pay anything until we can deliver your materials.

4.1 Custom Cloning

Fastest turn-around time of 1 ~ 2 weeks at \$350.00 (material cost), no matter what insert sizes are.

4.2 Site-direct Mutagenesis

Fastest turn-around time of 2 ~ 3 weeks at \$350.00 material cost, no matter how big insert sizes are.

4.3 Lentivirus

The most affordable lentivirus packaging system on earth starting at \$450.00/10 ml lentivirus at a concentration 10^6 TU/ml.

4.4 Adenovirus

Custom adenovirus for only \$650.00, cheaper than an adenovirus cloning kit on the market.

4.5 Stable Cell Line

The most reliable service for any stable cell line generation.

4.6 Custom Antibody

Both polyclonal and monoclonal antibody services are of great quality, beating any price quote from competitors.

4.7 Recombinant Protein Production

We are the leader in both bacteria and mammalian cell expression. Give us your project and we will deliver.

4.8 Custom Cell Immortalization

The leading service provider in the world for custom cell immortalization.

4.9 miRNA Profiling

Providing accurate service at material cost.

Product Name	Description	Size	Cat. No.	Price
Custom miRNA Profiling Service (Human)	1034 Mature Human miRNAs Based on miRBASE 16.0	Each Sample	C201	\$550.00
Custom miRNA Profiling Service (Mouse)	730 Mature Mouse miRNAs Based on miRBASE 16.0	Each Sample	C202	\$550.00
miRNA Profiling Kit (Human)	Human miRNAs Based on miRBASE 16.0	3 rxns	MA003	\$950.00
miRNA Profiling Kit (Mouse)	Mouse miRNAs Based on miRBASE 16.0	3 rxns	MA004	\$950.00

4.10 Whole Genome Sequencing

Standing confidently against the biggest names in the current market for demanding, high fidelity PCR applications and Whole Genome Sequencing, **abm**'s enzymes offer the highest accuracy rates available on the market with superior sensitivity and yield and extremely robust performance. To further harness the benefits of their proprietary enzyme technology, **abm** has established the Whole Genome Sequencing service below \$1000 not available from other providers.

4.11 Plasmid Sequencing

Accuracy is of paramount significance to any form of genetic materials used in life science research. Due to our large volume of sequencing, **abm** is able to offer quality plasmid sequencing at the most competitive price of \$3.50 per run and a flat rate of \$2.75 will be applicable for any order over 10 runs.

4.12 Oligo Synthesis

With the most advanced oligo synthesizer, **abm** is able to offer a deep, flat discount oligo synthesis at a unmatched \$0.15/bp in addition to other special oligo synthesis like fluorescent conjugated primers or Taqman probes.

Applied Biological Materials Inc.

Telephone:

(8:30am - 4:30pm PST, Mon - Fri)
Toll Free: 1-866-757-2414
Local: 604-247-2416
Fax: 604-247-2414 (24Hr.)

Email:

General Information: info@abmGood.com
Order: order@abmGood.com
Technical Support: technical@abmGood.com
Business Development: bd@abmGood.com

Address:

Suite #8-13520 Crestwood Place
Richmond, BC, Canada V6V 2G2

Website: www.abmGood.com

China

南京爱必梦应用生物材料有限公司

电话:

传真:

电邮: order@biocompete.com

网站: www.biocompete.com

北京艾可萃生物科技有限公司

电话: 64979359

传真: 64979359

电邮: tech_service@actbio.com.cn

网站: www.actbio.com.cn

重庆探生科技有限公司

电话: (23)68827630

400-602-3160

传真: (23)68104720

电邮: order@biom.cn

网站: www.biom.cn

Distributors

North America

Canada/USA

Applied Biological
Materials Inc.
Tel: 604-247-2416
order@abmGood.com

Mexico

Proveedor de Laboratorios S.A.C.V.
Tel: 52-33-38488910
prolab@pro-lab.com.mx

Puerto Rico

AVP Caribe
Tel: 787-892-0047
Fax: 787-264-3816
jcgonzalez@avpcaribe.com

Asia

South Korea

ITSBio, Inc.
Tel: 82-2-3462-8658
Fax: 82-2-3462-8659
bclim@itsbio.co.kr

Taiwan

Bio Pioneer Tech Co., Ltd.
Tel: 886-2-8660-9496
Fax: 886-2-8660-9342
tding.science@msa.hinet.net

Japan

Cosmo Bio Co. Ltd.
Tel: 03-5632-9610/9620
Fax: 03-5632-9619
mail@cosmobio.co.jp

India

Life Expression India
Tel: 91-998-698-8207
info@life-expression.com

Israel

BioConsult
Tel: 972-2-566-7043
Fax: 972-2-566-2790
sales@bioconsult.co.il

Singapore

Bio-REV PTE
Tel: 65-6273-3022
Fax: 65-6273-3020
allan@bio-rev.com

Europe

United Kingdom

NBS Biologicals Ltd.
Tel: 44-1480-433875
Fax: 44-1480-459868
info@nbsbio.co.uk

Belgium

Gentaur
Tel: 32-2-732-5688
Fax: 32-2-732-4414
ea@gentaur.com

France

Euromedex
Tel: 33-03-88-18-07-27
Fax: 33-03-88-18-07-28
info@euromedex.com

Germany

BioCat GmbH
Tel: 49-6221-714-1516
Fax: 49-6221-714-1529
info@biocat.com

Italy

Microtech s.r.l.
Tel: 39-0816-107435
Fax: 39-0816-107431
microtech@microtech.eu

Austria

THP Medical Products
Tel: 43-1-292-8280
Fax: 43-1-292-8280-88
office@thp.at