

## STAT1(Phospho-Tyr701) Antibody

Catalog No: #11044

Package Size: #11044-1 50ul #11044-2 100ul

## Description

|                       |                                                                                                                                                                                                                                                                                         |
|-----------------------|-----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Product Name          | STAT1(Phospho-Tyr701) Antibody                                                                                                                                                                                                                                                          |
| Host Species          | Rabbit                                                                                                                                                                                                                                                                                  |
| Clonality             | Polyclonal                                                                                                                                                                                                                                                                              |
| Purification          | Antibodies were produced by immunizing rabbits with synthetic phosphopeptide and KLH conjugates.<br>Antibodies were purified by affinity-chromatography using epitope-specific phosphopeptide. Non-phospho specific antibodies were removed by chromatography using non-phosphopeptide. |
| Applications          | WB IHC                                                                                                                                                                                                                                                                                  |
| Species Reactivity    | Hu Ms                                                                                                                                                                                                                                                                                   |
| Specificity           | The antibody detects endogenous level of STAT1 only when phosphorylated at tyrosine 701.                                                                                                                                                                                                |
| Immunogen Type        | Peptide-KLH                                                                                                                                                                                                                                                                             |
| Immunogen Description | Peptide sequence around phosphorylation site of tyrosine 701 (T-G-Y(p)-I-K) derived from Human STAT1.                                                                                                                                                                                   |
| Conjugates            | Unconjugated                                                                                                                                                                                                                                                                            |
| Target Name           | STAT1                                                                                                                                                                                                                                                                                   |
| Modification          | Phospho                                                                                                                                                                                                                                                                                 |
| Other Names           | CANDF7, ISGF-3, STAT91                                                                                                                                                                                                                                                                  |
| Accession No.         | Swiss-Prot: P42224NCBI Protein: NP_009330.1                                                                                                                                                                                                                                             |
| SDS-PAGE MW           | 84,91kd                                                                                                                                                                                                                                                                                 |
| Concentration         | 1.0mg/ml                                                                                                                                                                                                                                                                                |
| Formulation           | Supplied at 1.0mg/mL in phosphate buffered saline (without Mg <sup>2+</sup> and Ca <sup>2+</sup> ), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.                                                                                                                            |
| Storage               | Store at -20°C for long term preservation (recommended). Store at 4°C for short term use.                                                                                                                                                                                               |


## Application Details

Predicted MW: 84,91kd


Western blotting: 1:500~1:1000

Immunohistochemistry: 1:50~1:100


## Images



Immunohistochemical analysis of paraffin-embedded human breast carcinoma tissue using STAT1(Phospho-Tyr701) Antibody #11044(left) or the same antibody preincubated with blocking peptide(right).



Western blot analysis of extracts from MEF cells untreated or treated with interferon-I (IFN $\gamma$ ) using STAT1 (Phospho-Tyr701) Antibody #11044.



Western blot analysis of extracts from HeLa cells, treated with IFN $\alpha$  or calf intestinal phosphatase (CIP), using STAT1 (Phospho-Tyr701) Antibody #11044.

## Background

Signal transducer and activator of transcription that mediates signaling by interferons (IFNs). Following type I IFN (IFN- $\alpha$  and IFN- $\beta$ ) binding to cell surface receptors, Jak kinases (TYK2 and JAK1) are activated, leading to tyrosine phosphorylation of STAT1 and STAT2. The phosphorylated STATs dimerize, associate with ISGF3G/IRF-9 to form a complex termed ISGF3 transcription factor, that enters the nucleus. ISGF3 binds to the IFN stimulated response element (ISRE) to activate the transcription of interferon stimulated genes, which drive the cell in an antiviral state. In response to type II IFN (IFN- $\gamma$ ), STAT1 is tyrosine- and serine-phosphorylated. It then forms a homodimer termed IFN- $\gamma$ -activated factor (GAF), migrates into the nucleus and binds to the IFN  $\gamma$  activated sequence (GAS) to drive the expression of the target genes, inducing a cellular antiviral state.

Heim M H, (1999) J Recept Signal Transduct Res. 19: 75-120.

Durbin J E, et al. (1996) Cell. 84: 443-450.

Meraz M A, et al. (1996) Cell. 84: 431-442.

Wakao H, et al. (1994) EMBO J. 13: 2182-2191.

Demoulin J, B. et al. (1999) J Biol Chem. 274: 25855-258

Ihle J N, et al. (1994) Trends Biochem Sci. 19: 222-227.

## Published Papers

C.-C.E. Lan, C.-S. Wu, S.-M. Huang et al., High-glucose environment reduces human  $\beta$ -defensin-2 expression in human keratinocytes: implications for

poor diabetic wound healing., British Association of Dermatologists, 166(6):1221η— C1229(2012)

[PMID:22283836](#)

Hirokazu Hara, Yoko Nakamura, Masayuki Ninomiya et al., Inhibitory effects of chalcone glycosides isolated from *Brassica rapa* L.

and their synthetic derivatives on LPS-induced NO production in microglia., *Bioorganic & Medicinal Chemistry*, 19(18):5559-5568(2011)

[PMID:21856162](#)

Hsin-Chien Chen, Hsin-I Ma, Huey-Kang Sytwu et al., Neural Stem Cells Secrete Factors That Promote Auditory Cell Proliferation Via a Leukemia Inhibitory Factor Signaling Pathway, *Journal of Neuroscience Research*, 88:3308η— C3318(2010)

[PMID:20882565](#)

Jelke J. Fros, Wen Jun Liu, Natalie A. Prow et al., Chikungunya virus nonstructural protein 2 inhibits type I/II interferon-stimulated JAK-STAT signaling., *Journal of Virology*, 84(20):10877-10887(2010)

[PMID:20686047](#)

et al., Danzhi Jiangtang Capsule ameliorates kidney injury via inhibition of the JAK-STAT signaling pathway and increased antioxidant capacity in STZ-induced diabetic nephropathy rats. In *Biosci Trends*. On 2019 Jan 22 by Sun M, Bu W et al.. PMID: 30606979 , , (2018)

[PMID:30606979](#)

et al., Chronic Inflammation In Response to Injury: Retention of Myeloid Cells in Injured Tissue Is Driven by Myeloid Cell Intrinsic Factors. In *J Invest Dermatol* on 2019 Jan 28 by Torbica T, Wicks K, et al.. PMID:30703358, , (2019)

[PMID:30703358](#)

et al., A thiazole-derived oridonin analogue exhibits antitumor activity by directly and allosterically inhibiting STAT3. In *J Biol Chem* on 2019 Nov 15; by Shen X, Zhao L, et al.. PMID:31594861, , (2019)

[PMID:31594861](#)

et al., The C-terminal Domain of Chikungunya Virus nsP2 Independently Governs Viral RNA Replication, Cytopathicity, and Inhibition of Interferon Signaling. In *J Virol* on 2013 Sep by Jelke J Fros, Erika van der Maten, et al.. PMID:23864632, , (2013)

[PMID:23864632](#)

et al., High-glucose Environment Reduces Human  $\alpha$ -Y-defensin-2 Expression in Human Keratinocytes: Implications for Poor Diabetic Wound Healing. In

*Br J Dermatol* on 2012 Jun by C-C E Lan, C-S Wu, et al.. PMID:22283836, , (2012)

[PMID:22283836](#)

et al., Neural stem cells secrete factors that promote auditory cell proliferation via a leukemia inhibitory factor signaling pathway. In *J Neurosci Res* on 2010 Nov 15 by Hsin-Chien Chen, Hsin-I Ma, et al.. PMID: 20882565, , (2010)

[PMID:20882565](#)

et al., Chikungunya virus nonstructural protein 2 inhibits type I/II interferon-stimulated JAK-STAT signaling. In *J Virol* on 2010 Oct by Jelke J Fros, Wen Jun Liu, et al.. PMID:

20686047, , (2010)

[PMID:20686047](#)

et al., Inhibitory effects of chalcone glycosides isolated from *Brassica rapa* L. and their synthetic derivatives on LPS-induced NO production in microglia. In *Bioorg Med Chem* on 2011 Sep 15 by Hirokazu Hara, Yoko Nakamura, et al.. PMID:21856162, , (2011)

[PMID:21856162](#)

et al., DJC suppresses Advanced Glycation End Products-Induced JAK-STAT signaling and ROS in mesangial cells. In *Evid Based Complement Alternat Med* on 2017 May 29 by Min Sun, Yan Li, et al.. PMID: 28630633, , (2017)

[PMID:28630633](#)

et al., Inhibitory Effects of Chalcone Glycosides Isolated From *Brassica Rapa* L. 'Hidabeni' and Their Synthetic Derivatives on LPS-induced NO Production in Microglia. In *Bioorg Med Chem* on 2011 Sep 15 by Hirokazu Hara, Yoko Nakamura, et al.. PMID: 21856162, , (2017)

[PMID:21856162](#)

el at., Newly Synthesized 1,3-bis(4-iodophenyl)-5-hydroxy-4-methyl-2-pyrazolin-5-one Chalcone Derivatives Potently Suppress LPS-Induced NO Production via Inhibition of STAT1, but Not NF- $\kappa$ B, JNK, and p38, Pathways in Microglia. In Biol Pharm Bull On 2014 by Hirokazu Hara , Ryoko Ikeda et al..PMID: 24882415, , (2014)

[PMID:24882415](#)

el at., Molecular analysis of curcumin-induced polarization of murine RAW264.7 macrophages. In J Cardiovasc Pharmacol on 2014 Jun by Fangyuan Chen , Ning Guo et al..PMID:24709638, , (2014)

[PMID:24709638](#)

el at., Expression of interferon effector gene SART1 correlates with interferon treatment response against hepatitis B infection. In Mediators Inflamm on 2016 by Yong Li, Chuanlong Zhu et al..PMID: 28077916, , (2016)

[PMID:28077916](#)

el at., Mycobacterium tuberculosis EspB protein suppresses interferon- $\gamma$ -induced autophagy in murine macrophages. In J Microbiol Immunol Infect on 2016 Dec by DanDan Huang, Lang Bao et al..PMID:25641595, , (2016)

[PMID:25641595](#)

el at., Costunolide improved dextran sulfate sodium-induced acute ulcerative colitis in mice through NF- $\kappa$ B, STAT1/3, and Akt signaling pathways. In Int Immunopharmacol on 2020 Jul by Fan Xie, Hai Zhang, et al..PMID:32413737, , (2020)

[PMID:32413737](#)

Ighodaro Igbe;XiaoFei Shen;Wei Jiao;Zhe Qiang;Teng Deng;Sheng Li;WanLi Liu;HanWei Liu;GuoLin Zhang;Fei Wang el at., Dietary quercetin potentiates the antiproliferative effect of interferon- $\alpha$  in hepatocellular carcinoma cells through activation of JAK/STAT pathway signaling by inhibition of SHP2 phosphatase, , (2017)

[PMID:29371942](#)

Rui Yang;Meng-Yan Deng;Lu-Kui Yang;Guan-Zhe Wang;Jun Ma;Qiang Wen;Na Gao;Hai-Ling Qiao el at., Identification of Cytochrome P450 2E1 as a Novel Target in Neuroinflammation and Development of Its Inhibitor Q11 as a Treatment Strategy., , (2025)

[PMID:40122152](#)

---

Note: This product is for in vitro research use only and is not intended for use in humans or animals.