

FoxO3a (Phospho-Ser574) Antibody

Catalog No: #12874

Package Size: #12874-1 50ul #12874-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	FoxO3a (Phospho-Ser574) Antibody
Host Species	Rabbit
Clonality	Polyclonal
Applications	WB
Species Reactivity	Hu Ms Rt
Specificity	Phospho-FoxO3a (S574) Antibody detects endogenous levels of FoxO3a only when phosphorylated at S574
Immunogen Type	Peptide-KLH
Immunogen Description	A synthesized peptide derived from human FoxO3a (Phospho-Ser574)
Conjugates	Unconjugated
Other Names	AF6q21 antibody AF6q21 protein antibody DKFZp781A0677 antibody FKHR2 antibody FKHRL 1 antibody FKHRL1 antibody FKHRL1P2 antibody Forkhead (Drosophila) homolog (rhabdomyosarcoma) like 1 antibody Forkhead box O3 antibody Forkhead box O3A antibody Forkhead box protein O3 antibody Forkhead box protein O3A antibody Forkhead Drosophila homolog of in rhabdomyosarcoma like 1 antibody Forkhead homolog (rhabdomyosarcoma) like 1 antibody Forkhead in rhabdomyosarcoma like 1 antibody Forkhead in rhabdomyosarcoma-like 1 antibody FOX O3A antibody FOXO2 antibody foxo3 antibody FOXO3_HUMAN antibody FOXO3A antibody MGC12739 antibody MGC31925 antibody
Accession No.	Swiss-Prot#:O43524 NCBI Gene ID2309
Calculated MW	97
Concentration	1.0mg mL
Formulation	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺) pH 7.4 150mM NaCl 0.02% sodium azide and 50% glycerol.

Storage

Store at -20°C

Application Details

WB dilution:1:1000

Published Papers

Cui Ran-Ji;Gao Shuo-Hui;Guo Qiu-Shi;Li Bing-Jin;Sun Li-Rui;Sun Zhi-Hui;Yang Wei;Zhang Hong-Mei;Zhou Wei el at., Modulation of Multiple Signaling Pathways of the Plant-Derived Natural Products in Cancer, , (2019)

PMID:31781485

Note: This product is for in vitro research use only and is not intended for use in humans or animals.