

MOB1A (Phospho-Tyr26) Antibody

Catalog No: #12878

Package Size: #12878-1 50ul #12878-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	MOB1A (Phospho-Tyr26) Antibody
Host Species	Rabbit
Clonality	Polyclonal
Applications	WB
Species Reactivity	Hu Ms Rt
Specificity	Phospho-MOB1A (Y26) Antibody detects endogenous levels of MOB1A only when phosphorylated at Y26
Immunogen Type	Peptide-KLH
Immunogen Description	A synthesized peptide derived from human MOB1A (Phospho-Tyr26)
Conjugates	Unconjugated
Other Names	MATS 2 antibody MATS2 antibody MGC33910 antibody Mob 1A antibody Mob 1B antibody MOB 4A antibody MOB kinase activator 1B antibody Mob1 homolog 1A antibody MOB1 Mps One Binder homolog B antibody MOB1 Mps one binder kinase activator like 1A antibody Mob1A antibody Mob1B antibody MOBKL 1A antibody MOBKL1A antibody MOL1A_HUMAN antibody Mps one binder kinase activator like 1A antibody Mps one binder kinase activator-like 1A antibody Protein Mob4A antibody
Accession No.	Swiss-Prot#:Q7L9L4 NCBI Gene ID92597
Calculated MW	25
Concentration	1.0mg mL
Formulation	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺) pH 7.4 150mM NaCl 0.02% sodium azide and 50% glycerol.
Storage	Store at -20°C

Application Details

WB dilution:1:1000

Published Papers

el et., A Platform of Synthetic Lethal Gene Interaction Networks Reveals that the GNAQ Uveal Melanoma Oncogene Controls the Hippo Pathway through FAK. In Cancer Cell on 2019 Mar 18 by Feng X, Arang N, et al.. PMID:30773340, , (2019)

[PMID:30773340](#)

Note: This product is for in vitro research use only and is not intended for use in humans or animals.