

Vimentin (Phospho-S39) Rabbit mAb

Catalog No: #14301

Package Size: #14301-1 50ul #14301-2 100ul


Description

Product Name	Vimentin (Phospho-S39) Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Clone No.	SR1953
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB
Species Reactivity	Human Mouse Rat
Specificity	Phospho-Vimentin (S39) Antibody detects endogenous levels of total Phospho-Vimentin (S39)
Immunogen Description	A synthesized peptide derived from human Phospho-Vimentin (S39)
Conjugates	Unconjugated
Other Names	VIM; VIME; Vimentin;
Accession No.	P08670
Calculated MW	Predicted band size: 54 kDa
SDS-PAGE MW	Observed band size: 54 kDa
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at -20°C for long term preservation (recommended). Store at 4°C for short term use.

Application Details

WB: 1:500-1:2000

Images



All lanes : Vimentin (Phospho-S39) Rabbit mAb at 1/1k dilution
Lane 1 : HeLa whole cell lysates
Lane 2 : HeLa treated with 100nM Calyculin A for 30min whole cell lysates
Lysates/proteins at 20 µg per lane.
Secondary All lanes : Goat Anti-Rabbit IgG H&L (HRP) at 1/20000 dilution
Predicted band size: 54 kDa
Observed band size: 54 kDa
Exposure time: 2 seconds

Background

Phospho-Vimentin (S39) Antibody detects endogenous levels of total Phospho-Vimentin (S39)

Note: This product is for in vitro research use only and is not intended for use in humans or animals.