

Integrin b3(Ab-785) Antibody

Catalog No: #21274

Package Size: #21274-1 50ul #21274-2 100ul

Description


Product Name	Integrin b3(Ab-785) Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Antibodies were produced by immunizing rabbits with synthetic peptide and KLH conjugates. Antibodies were purified by affinity-chromatography using epitope-specific peptide.
Applications	WB
Species Reactivity	Hu Ms Rt
Specificity	The antibody detects endogenous level of total Integrin b3 protein.
Immunogen Type	Peptide-KLH
Immunogen Description	Peptide sequence around aa.783~787 (I-T-Y-R-G) derived from Human Integrin B3/ITGB3.
Conjugates	Unconjugated
Target Name	Integrin b3
Other Names	CD61 antigen; GP3A; GPIIIa; ITB3; Platelet membrane glycoprotein IIIa
Accession No.	Swiss-Prot: P05106NCBI Protein: NP_000203.2
Concentration	1.0mg/ml
Formulation	Supplied at 1.0mg/mL in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at -20°C for long term preservation (recommended). Store at 4°C for short term use.

Application Details

Predicted MW: 110kd

Western blotting: 1:500~1:1000

Images



Western blot analysis of extracts from HL60 and HT29 cells using Integrin b3(Ab-785) Antibody #21274.

Background

The Itg-Beta3 gene encodes the beta subunit of the platelet membrane adhesive protein receptor complex Itg-Beta3 protein, which belongs to the Integrin class of cell adhesion molecule receptors that share a common heterotrimeric structure with α and Beta subunits. It is a common Beta subunit of the platelet complex and of two other Integrins (Fibronectin and Vitronectin Receptors), which have distinct α subunits. Itg-Beta3 has been implicated in a wide variety of functions, including platelet aggregation and thrombosis and implantation, placentation, angiogenesis, bone remodeling, and tumor progression. Glanzmann Thrombasthenia can result from defects in the genes for either the Itg- α 2B or the Itg-Beta3 subunit
Sujoy Bhattacharya, et al. (2006) Biochem J. August 1; 397(Pt 3): 437

Published Papers

Ruifei Wang, Jiajia Bi, Khamal Kwesi Ampah et al., Lipid raft regulates the initial spreading of melanoma A375 cells by modulating α 1 integrin clustering., The International Journal of Biochemistry & Cell Biology., 45(8):1679-1689(2013)

[PMID:23665237](#)

et al., Lipid raft regulates the initial spreading of melanoma A375 cells by modulating α 1Y1 integrin clustering. In Int J Biochem Cell Biol on 2013 Aug by Ruifei Wang, Jiajia Bi, et al.. PMID: 23665237, , (2013)

[PMID:23665237](#)

Note: This product is for in vitro research use only and is not intended for use in humans or animals.