

LKB1 Antibody


Catalog No: #24473

Package Size: #24473 100ul

Description

Product Name	LKB1 Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Affinity chromatography purified via peptide column
Applications	ELISA WB
Species Reactivity	Hu Ms
Immunogen Type	Peptide
Immunogen Description	Raised against a 15 amino acid peptide from near the carboxy terminus of human LKB1.
Conjugates	Unconjugated
Target Name	LKB1
Other Names	Serine-threonine protein kinase 11, STK11, NT-REN-19 antigen
Accession No.	Swiss-Prot:Q15831Gene ID:6794
Concentration	1mg/ml
Formulation	Supplied in PBS containing 0.02% sodium azide.
Storage	Can be stored at -20°C, stable for one year. As with all antibodies care should be taken to avoid repeated freeze thaw cycles. Antibodies should not be exposed to prolonged high temperatures.

Images



Western blot analysis of LKB1 in PC-3 cell lysate with LKB1 antibody at 1 ug/mL in the (A) absence or (B) presence of blocking peptide.

Background

The LKB1 serine/threonine protein kinase was initially identified as a tumor suppressor gene mutated in human Peutz-Jeghers syndrome (PJS), a condition resulting in the growth of numerous intestinal polyps classed as hamartomas. LKB1 exists as a heterotrimeric complex with two other proteins, Ste20-related adaptor protein (STRAD) and MO25. Together, this complex can phosphorylate and activate the AMP-activate protein kinase (AMPK). Following AMPK activation by LKB1, AMPK then phosphorylates TSC1 and TSC2, key components of the metabolism-regulating TOR signaling pathway, which antagonizes the activation for the TOR pathway. LKB1 has also been shown to play a fundamental role in controlling the spatial orientation of structures required to maintain an ordered, polarized epithelium.

Published Papers

el at., Artemisinin conferred cytoprotection to human retinal pigment epithelial cells exposed to amiodarone-induced oxidative insult by activating the CaMKK2/AMPK/Nrf2. In J Transl Med on 2024 Sep 16 by Chao Yang, Xia Zhao,et al..PMID:39285426, , (2024)

[PMID:39285426](#)

Note: This product is for in vitro research use only and is not intended for use in humans or animals.