

PerCP-Cy5.5 anti-mouse CD185

Catalog No: #26844

Package Size: #26844 100Tests



Description

Product Name	PerCP-Cy5.5 anti-mouse CD185
Host Species	Rat
Clonality	Monoclonal
Isotype	Rat IgG2b, k
Applications	FC
Species Reactivity	Mouse
Conjugates	Unconjugated
Storage	Store at 4°C in dark

Application Details

FC 5μl/time

Note: This product is for in vitro research use only and is not intended for use in humans or animals.