

PE anti-mouse F4/80

Catalog No: #26848

Package Size: #26848 100 Tests



SAB[®]
Signalway Antibody

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com

Description

Product Name	PE anti-mouse F4/80
Isotype	rat IgG2b
Applications	FC
Species Reactivity	Mouse
Conjugates	Unconjugated
Other Names	EMR1, Ly71
Formulation	0.2%BSAε 0.1%Proclin 300 Phosphate buffered solutionoO PH7.2-7.4
Storage	Store at 4°C in dark

Application Details

FC 5μl/time

Note: This product is for in vitro research use only and is not intended for use in humans or animals.