

Alexa Fluor 700 anti-mouse CD86

Catalog No: #26849

Package Size: #26849 100 Tests



Description

|                    |                                                                  |
|--------------------|------------------------------------------------------------------|
| Product Name       | Alexa Fluor 700 anti-mouse CD86                                  |
| Isotype            | Rat IgG2a, k                                                     |
| Applications       | FC                                                               |
| Species Reactivity | Mouse                                                            |
| Conjugates         | Unconjugated                                                     |
| Other Names        | B7-2, B70, Ly-58                                                 |
| Formulation        | 0.2%BSAε 0.1%Proclin 300 Phosphate buffered solutionoO PH7.2-7.4 |
| Storage            | Store at 4°C in dark                                             |

Application Details

|                |
|----------------|
| FC 5&mu;l/time |
|----------------|

Note: This product is for in vitro research use only and is not intended for use in humans or animals.