

APC anti-mouse CD206

Catalog No: #26850

Package Size: #26850 100 Tests



SAB[®]
Signalway Antibody

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com

Description

Product Name	APC anti-mouse CD206
Applications	FC
Species Reactivity	Mouse
Conjugates	Unconjugated
Formulation	0.2%BSAε 0.1%Proclin 300 Phosphate buffered solutionoO PH7.2-7.4
Storage	Store at 4°C in dark

Application Details

FC 5μl/time

Note: This product is for in vitro research use only and is not intended for use in humans or animals.