

## TGFβ1 Antibody

Catalog No: #39303

Package Size: #39303 100ul


## Description

|                       |                                                                                                                                                                                                                                                                                 |
|-----------------------|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Product Name          | TGFβ1 Antibody                                                                                                                                                                                                                                                                  |
| Host Species          | Rabbit                                                                                                                                                                                                                                                                          |
| Clonality             | Polyclonal                                                                                                                                                                                                                                                                      |
| Purification          | purified                                                                                                                                                                                                                                                                        |
| Applications          | ELISA WB IHC                                                                                                                                                                                                                                                                    |
| Species Reactivity    | Hu Ms Rt                                                                                                                                                                                                                                                                        |
| Specificity           | Not yet tested in other species.                                                                                                                                                                                                                                                |
| Immunogen Type        | Peptide-KLH                                                                                                                                                                                                                                                                     |
| Immunogen Description | Synthetic peptide conjugated to KLH                                                                                                                                                                                                                                             |
| Conjugates            | Unconjugated                                                                                                                                                                                                                                                                    |
| Target Name           | TGFβ1                                                                                                                                                                                                                                                                           |
| Concentration         | 1.0mg/ml                                                                                                                                                                                                                                                                        |
| Formulation           | PBS, pH 7.4 with 0.02% Sodium Azide                                                                                                                                                                                                                                             |
| Storage               | This product is stable for several weeks at 4°C as an undiluted liquid. Dilute only prior to immediate use. For extended storage, aliquot contents and freeze at -20°C or below. Avoid cycles of freezing and thawing.<br>Expiration date is one (1) year from date of receipt. |


## Application Details

ELISA WB IHC

## Images



Immunohistochemical staining of formalin-fixed paraffin-embedded heart showing cytoplasmic staining with anti-TGF beta1 antibody dilution:1:100



## Published Papers

JingBei Zhang;HongLiu Jin;XiaoYing Feng;Senling Feng;WenTing Zhu;HongMei Nan;ZhongWen Yuan et al., The combination of Lonicerae Japonicae Flos and Forsythiae Fructus herb-pair alleviated inflammation in liver fibrosis, , (2022)

[PMID:36059967](#)

Note: This product is for in vitro research use only and is not intended for use in humans or animals.