

LIN41 Antibody

Catalog No: #39754

Package Size: #39754 100ul


Description

Product Name	LIN41 Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Affinity purified
Applications	ELISA WB
Specificity	Human and predicted against mouse and rat.
Immunogen Type	Peptide-KLH
Immunogen Description	Rabbit polyclonal LIN41 (1) antibody was raised against a peptide from N terminal residues of human LIN41 protein.
Conjugates	Unconjugated
Target Name	LIN41
Concentration	1.0mg/ml
Formulation	PBS, pH 7.4 with 0.02% Sodium Azide
Storage	This product is stable for several weeks at 4°C as an undiluted liquid. Dilute only prior to immediate use. For extended storage, aliquot contents and freeze at -20°C or below. Avoid cycles of freezing and thawing. Expiration date is one (1) year from date of receipt.

Application Details

ELISA: 1:20000-1:80000 WB: 1:1000-1:2000.

Images



Published Papers

el at., Ilex pubescens inhibits pyroptosis post-myocardial infarction through suppression of the ROS/NLRP3 pathway, , (2024)

PMID:

Note: This product is for in vitro research use only and is not intended for use in humans or animals.