

GGTLC1 Antibody

Catalog No: #40020

Package Size: #40020 100ul



SAB[®]
Signalway Antibody

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com


Description

Product Name	GGTLC1 Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Affinity purified
Applications	ELISA WB IHC
Species Reactivity	Hu Rt Ms
Immunogen Type	Recombinant Protein
Immunogen Description	Rabbit polyclonal GGTLC1(1) antibody was raised against a recombinant human GGTLC1 protein 108-221aa(BC040904).
Conjugates	Unconjugated
Target Name	GGTLC1
Concentration	1.0mg/ml
Formulation	PBS, pH 7.4 with 0.02% sodium azide.
Storage	This product is stable for several weeks at 4°C as an undiluted liquid. Dilute only prior to immediate use. For extended storage, aliquot contents and freeze at -20°C or below. Avoid cycles of freezing and thawing. Expiration date is one (1) year from date of receipt.

Application Details

ELISA: 1:20000-1:80000 WB: 1:200-1:1000 IHC: 1:100-1:500

Images



Published Papers

el at., GGTLC1 knockdown inhibits the progression of endometrial cancer by regulating the TGF-β/Smad signaling pathway. In Heliyon on 2024 May 25 by Xiaolin Tang, Xuehan Bi,et al..PMID:338841516, , (2024)

[PMID:338841516](https://pubmed.ncbi.nlm.nih.gov/338841516/)

Note: This product is for in vitro research use only and is not intended for use in humans or animals.