

ERMP1 Antibody

Catalog No: #40050

Package Size: #40050 100ul


Description

Product Name	ERMP1 Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Affinity purified
Applications	ELISA WB
Species Reactivity	Hu Ms Rt
Immunogen Type	Recombinant Protein
Immunogen Description	Rabbit polyclonal ERMP1(1)antibody was raised against a recombinant humanERMP1 protein 204-365 aa aa(BC031630)
Conjugates	Unconjugated
Target Name	ERMP1
Concentration	2.8mg/ml
Formulation	PBS, pH 7.4 with 0.02% sodium azide.
Storage	This product is stable for several weeks at 4°C as an undiluted liquid. Dilute only prior to immediate use. For extended storage, aliquot contents and freeze at -20°C or below. Avoid cycles of freezing and thawing. Expiration date is one (1) year from date of receipt.

Application Details

ELISA: 1:20000-1:80000 WB: 1:200-1:1000

Images



Published Papers

el at., miR-9-5p attenuates ischemic stroke through targeting ERMP1-mediated endoplasmic reticulum stress. In Acta Histochem on 2019 Nov by Chi L, Jiao D,et al..PMID:31500865, , (2019)

[PMID:31500865](https://pubmed.ncbi.nlm.nih.gov/31500865/)

Note: This product is for in vitro research use only and is not intended for use in humans or animals.