

IFNA7 Antibody

Catalog No: #45623

Package Size: #45623-1 50ul #45623-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	IFNA7 Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Immunogen affinity purified
Applications	WB
Species Reactivity	Hu
Specificity	IFNA7 Antibody detects endogenous levels of total IFNA7
Immunogen Type	Peptide
Immunogen Description	A synthesized peptide derived from human IFNA7
Conjugates	Unconjugated
Target Name	IFNA7
Other Names	IFNA7, IFN-alpha 7, IFN-alpha-7, IFN-alpha-J1, IFN-alphaJ, IFNA-J, Interferon alpha-J, Interferon alpha-J1, Interferon alpha-7, LeIF J, Interferon, alpha 7
Accession No.	Swiss-Prot#: P01567NCBI Gene ID:3444
Calculated MW	22 kDa
Concentration	1.0mg/mL
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at -20°C

Application Details

Western blotting: 1:1000-3000

Images



Western blot analysis IFNA7 using HepG2 whole cell lysates

Note: This product is for in vitro research use only and is not intended for use in humans or animals.