

c-Myc (Phospho-T58+S62) Antibody

Catalog No: #48685

Package Size: #48685-1 50ul #48685-2 100ul



Signalway Antibody

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com


Description

Product Name	c-Myc (Phospho-T58+S62) Antibody
Host Species	Recombinant Rabbit
Clonality	Monoclonal
Clone No.	22
Purification	ProA affinity purified
Applications	WB, ICC/IF, IP
Species Reactivity	Hu, Rt
Immunogen Description	Synthetic phospho-peptide corresponding to residues surrounding Thr58 and Ser62 of human c-Myc.
Conjugates	Unconjugated
Other Names	Avian myelocytomatosis viral oncogene homolog antibody bHLHe39 antibody c Myc antibody Class E basic helix-loop-helix protein 39 antibody MRTL antibody Myc antibody Myc protein antibody Myc proto oncogene protein antibody Myc proto-oncogene protein antibody myc related translation/localization regulatory factor antibody MYC_HUMAN antibody Myc2 antibody MYCC antibody Niard antibody Nird antibody Proto-oncogene c-Myc antibody Transcription factor p64 antibody v myc avian myelocytomatosis viral oncogene homolog antibody v myc myelocytomatosis viral oncogene homolog antibody
Accession No.	Swiss-Prot#:P01106
Calculated MW	49 kDa
Formulation	1*TBS (pH7.4), 1%BSA, 40%Glycerol. Preservative: 0.05% Sodium Azide.
Storage	Store at -20Λ C

Application Details

WB: 1:1,000 ICC: 1:50-1:200

Images



Western blot analysis of Phospho-c-Myc(T58+S62) on K562 cells lysates using anti-Phospho-c-Myc(T58+S62) antibody at 1/1,000 dilution.

Product Description

c-Myc (Phospho-T58+S62) Antibody

Note: This product is for in vitro research use only and is not intended for use in humans or animals.