

Ezrin Rabbit mAb

Catalog No: #49438

Package Size: #49438-1 50ul #49438-2 100ul



SAB[®]
Signalway Antibody

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com


Description

Product Name	Ezrin Rabbit mAb
Host Species	Recombinant Rabbit
Clonality	Monoclonal
Clone No.	JM10-65
Purification	ProA affinity purified
Applications	WB, IHC
Species Reactivity	Hu
Immunogen Description	recombinant protein
Conjugates	Unconjugated
Other Names	Villin 2 ezrin antibody CVIL antibody CVL antibody Cytovillin 2 antibody Cytovillin antibody DKFZp762H157 antibody Epididymis secretory protein Li 105 antibody EZR antibody EZRI_HUMAN antibody Ezrin antibody FLJ26216 antibody HEL S 105 antibody MGC1584 antibody p81 antibody VIL 2 antibody VIL2 antibody Villin 2 (ezrin) antibody Villin 2 antibody Villin-2 antibody Villin2 antibody
Accession No.	Swiss-Prot#:P15311
Calculated MW	80 kDa
Formulation	1*TBS (pH7.4), 1%BSA, 40%Glycerol. Preservative: 0.05% Sodium Azide.
Storage	Store at -20°C


Application Details

WB: 1:500-1:1000IHC: 1:50-1:200


Images




Western blot analysis of Ezrin on different cell lysates using anti-Ezrin antibody at 1/1,000 dilution. Positive control: Lane 1:Hela Lane 2: Jurkat Lane 3:A431




Immunohistochemical analysis of paraffin-embedded human lung tissue using anti-Ezrin antibody. Counter stained with hematoxylin.




Immunohistochemical analysis of paraffin-embedded human lung cancer tissue using anti-Ezrin antibody. Counter stained with hematoxylin.




Immunohistochemical analysis of paraffin-embedded human liver cancer tissue using anti-Ezrin antibody. Counter stained with hematoxylin.



Immunohistochemical analysis of paraffin-embedded human colon cancer tissue using anti-Ezrin antibody. Counter stained with hematoxylin.



Immunohistochemical analysis of paraffin-embedded human breast carcinoma tissue using anti-Ezrin antibody. Counter stained with hematoxylin.



Immunohistochemical analysis of paraffin-embedded human kidney tissue using anti-Ezrin antibody. Counter stained with hematoxylin.

References

Note: This product is for in vitro research use only and is not intended for use in humans or animals.