

Glutathione S Transferase mu 1 Rabbit mAb

Catalog No: #52314

Package Size: #52314-1 50ul #52314-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	Glutathione S Transferase mu 1 Rabbit mAb
Clonality	Monoclonal
Clone No.	S03-8G2
Isotype	Rabbit IgG
Purification	Affinity Purified
Applications	WB
Species Reactivity	Human,Mouse,Rat
Immunogen Description	A synthetic peptide of human Glutathione S Transferase mu
Conjugates	Unconjugated
Modification	Unmodification
Other Names	MU; H-B; GST1; GTH4; GTM1; MU-1; GSTM1-1; GSTM1a-1a; GSTM1b-1b
Accession No.	Swiss-Prot:P09488GenelD:2944
Calculated MW	Calculated MW: 26 kDa; Observed MW: 26 kDa
Concentration	0.3 mg/ml
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.

Application Details

WB 1:500-1:2000

Images



Western blot analysis of GSTM1 expression in SH-SY5Y cell lysate.

Background

Conjugation of reduced glutathione to a wide number of exogenous and endogenous hydrophobic electrophiles.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.