

NR1D1 Rabbit mAb

Catalog No: #52470

Package Size: #52470-1 50ul #52470-2 100ul


Description

|                       |                                                                                          |
|-----------------------|------------------------------------------------------------------------------------------|
| Product Name          | NR1D1 Rabbit mAb                                                                         |
| Host Species          | Recombinant Rabbit                                                                       |
| Clonality             | Monoclonal                                                                               |
| Clone No.             | S05-8B1                                                                                  |
| Isotype               | Rabbit IgG                                                                               |
| Purification          | Affinity Purified                                                                        |
| Applications          | WB IF                                                                                    |
| Species Reactivity    | Human,Mouse,Rat                                                                          |
| Immunogen Description | A synthetic peptide of human NR1D1                                                       |
| Conjugates            | Unconjugated                                                                             |
| Modification          | Unmodification                                                                           |
| Other Names           | EAR1; hRev; THRA1; THRAL; ear-1; REVERBA; REVERBalpha                                    |
| Accession No.         | Swiss-Prot:P20393GenelD:9572                                                             |
| Calculated MW         | Calculated MW: 67 kDa; Observed MW: 67 kDa                                               |
| Concentration         | 0.3 mg/ml                                                                                |
| Formulation           | 50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% BSA    |
| Storage               | Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles. |

Application Details


WB: 1/1000; ICC/IF: 1/20;

Images



Western blot detection of NR1D1 in K562,C6,3T3,Hela cell lysates using NR1D1 Rabbit mAb(1:1000 diluted).Predicted band size:67kDa.Observed band size:67kDa.

Immunofluorescence of NR1D1(green) in Hela cells using NR1D1 Rabbit mAb at dilution 1/50, and DAPI(blue)



## Background

Swiss-Prot Acc.P20393. Transcriptional repressor which coordinates circadian rhythm and metabolic pathways in a heme-dependent manner. Integral component of the complex transcription machinery that governs circadian rhythmicity and forms a critical negative limb of the circadian clock by directly repressing the expression of core clock components ARTNL/BMAL1, CLOCK and CRY1. Also regulates genes involved in metabolic functions, including lipid and bile acid metabolism, adipogenesis, gluconeogenesis and the macrophage inflammatory response. Acts as a receptor for heme which stimulates its interaction with the NCOR1/HDAC3 corepressor complex, enhancing transcriptional repression. Recognizes two classes of DNA response elements within the promoter of its target genes and can bind to DNA as either monomers or homodimers, depending on the nature of the response element. Binds as a monomer to a response element composed of the consensus half-site motif 5'-[A/G]GGTCA-3' preceded by an A/T-rich sequence (RevRE), or as a homodimer to a direct repeat of the core motif spaced by two nucleotides (RevDR-2). Acts as a potent competitive repressor of ROR alpha (RORA) function and regulates the levels of its ligand heme by repressing the expression of PPARGC1A, a potent inducer of heme synthesis. Regulates lipid metabolism by repressing the expression of APOC3 and by influencing the activity of sterol response element binding proteins (SREBPs); represses INSIG2 which interferes with the proteolytic activation of SREBPs which in turn govern the rhythmic expression of enzymes with key functions in sterol and fatty acid synthesis. Regulates gluconeogenesis via repression of G6PC and PEPCK and adipocyte differentiation via repression of PPARG. Regulates glucagon release in pancreatic alpha-cells via the AMPK-NAMPT-SIRT1 pathway and the proliferation, glucose-induced insulin secretion and expression of key lipogenic genes in pancreatic-beta cells. Positively regulates bile acid synthesis by increasing hepatic expression of CYP7A1 via repression of NR0B2 and NFIL3 which are negative regulators of CYP7A1. Modulates skeletal muscle oxidative capacity by regulating mitochondrial biogenesis and autophagy; controls mitochondrial biogenesis and respiration by interfering with the STK11-PRKAA1/2-SIRT1-PPARGC1A signaling pathway. Represses the expression of SERPINE1/PAI1, an important modulator of cardiovascular disease and the expression of inflammatory cytokines and chemokines in macrophages. Represses gene expression at a distance in macrophages by inhibiting the transcription of enhancer-derived RNAs (eRNAs). Plays a role in the circadian regulation of body temperature and negatively regulates thermogenic transcriptional programs in brown adipose tissue (BAT); imposes a circadian oscillation in BAT activity, increasing body temperature when awake and depressing thermogenesis during sleep. In concert with NR2E3, regulates transcriptional networks critical for photoreceptor development and function. In addition to its activity as a repressor, can also act as a transcriptional activator. In the ovarian granulosa cells acts as a transcriptional activator of STAR which plays a role in steroid biosynthesis. In collaboration with SP1, activates GJA1 transcription in a heme-independent manner.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.