

HIF-1 alpha Rabbit mAb

Catalog No: #56033

Package Size: #56033-1 50ul #56033-2 100ul



SAB[®]
Signalway Antibody

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com


Description

Product Name	HIF-1 alpha Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB IHC ICC/IF IP FC
Species Reactivity	Human
Specificity	HIF-1 alpha Antibody detects endogenous levels of total HIF-1 alpha
Immunogen Description	A synthesized peptide derived from human HIF-1 alpha
Conjugates	Unconjugated
Other Names	HIF1; MOP1; PASD8; bHLHe78; HIF-1alpha; HIF1-ALPHA; HIF1A
Accession No.	Uniprot:Q16665
Calculated MW	120kDa
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at +4ℳ C short term. Store at -20ℳ C long term. Avoid freeze / thaw cycle.


Application Details

WB:1:500~1:2000IHC:1:50~1:200ICC/IF:1:50~1:200IP:1:30FC:1:30


Images



Western blot analysis of HIF-1 alpha expression in Ramos cell lysate.



Immunohistochemical analysis of paraffin-embedded human kidney, using HIF-1 alpha Antibody.



Immunofluorescent analysis of HeLa cells, using HIF-1 alpha Antibody.

Published Papers

Jingyao Li;Huixi Yi;Yuanyuan Fu;Jiani Zhuang;Zhixiong Zhan;Liyu Guo;Ji Zheng;Xiyong Yu;Dong-Yang Zhang et al., Biodegradable iridium coordinated nanodrugs potentiate photodynamic therapy and immunotherapy of lung cancer., (2025)

[PMID:39488900](#)

Note: This product is for in vitro research use only and is not intended for use in humans or animals.