

CD3 zeta Rabbit mAb

Catalog No: #56124

Package Size: #56124-1 50ul #56124-2 100ul

Description


Product Name	CD3 zeta Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Clone No.	SR2324
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	ICC/IF FC
Species Reactivity	Human
Specificity	CD3 zeta Antibody detects endogenous levels of total CD3 zeta
Immunogen Description	A synthesized peptide derived from human CD3 zeta
Conjugates	Unconjugated
Other Names	CD3-zeta; T-cell receptor T3 zeta chain; T3Z; TCRZ; CD247;
Accession No.	Uniprot:P20963
Calculated MW	kDa
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at -20°C for long term preservation (recommended). Store at 4°C for short term use.

Application Details

ICC/IF:1:100~1:200

FC:1:20

Images



All lanes use the Antibody at 1:1K dilution for 1 hour at room temperature.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.