

VE Cadherin Rabbit mAb

Catalog No: #56125

Package Size: #56125-1 50ul #56125-2 100ul


Description

Product Name	VE Cadherin Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB ICC/IF IP
Species Reactivity	Mouse
Specificity	VE Cadherin Antibody detects endogenous levels of total VE Cadherin
Immunogen Description	A synthesized peptide derived from mouse VE Cadherin
Conjugates	Unconjugated
Other Names	Vascular endothelial cadherin; VE-cadherin; 7B4 antigen; CD144 antigen; CDH5;
Accession No.	Uniprot:P55284(mouse)
Calculated MW	90,125kDa
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at +4 Λ C short term. Store at -20 Λ C long term. Avoid freeze / thaw cycle.

Application Details

WB:1:500~1:1000ICC/IF:1:500~1:1000IP:1:50

Images



Western blot analysis of VE Cadherin expression in Mouse lung lysate.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.