

PGC1 alpha/beta Rabbit mAb

Catalog No: #56129

Package Size: #56129-1 50ul #56129-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	PGC1 alpha/beta Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB
Species Reactivity	Human Mouse Rat
Specificity	PGC1 alpha/beta Antibody detects endogenous levels of total PGC1 alpha/beta
Immunogen Description	A synthesized peptide derived from human PGC1 alpha/beta
Conjugates	Unconjugated
Other Names	PPARGC1A; LEM6; PGC1; PGC1A; PGC-1v; PPARGC1; PGC-1(alpha);
Accession No.	Uniprot:Q9UBK2/Q86YN6
Calculated MW	113kDa
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at +4ℳ C short term. Store at -20ℳ C long term. Avoid freeze / thaw cycle.

Application Details

WB:1:500~1:2000

Images



Western blot analysis of PGC1 alpha/beta expression in (1) HeLa cell lysate; (2) NIH 3T3 cell lysate; (3) C6 cell lysate.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.