

Substance P Rabbit mAb

Catalog No: #56143

Package Size: #56143-1 50ul #56143-2 100ul


Description

Product Name	Substance P Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	IHC
Species Reactivity	Human
Specificity	Substance P Antibody detects endogenous levels of total Substance P
Immunogen Description	A synthesized peptide derived from human Substance P
Conjugates	Unconjugated
Other Names	Protachykinin-1; PPT; Neurokinin A; Neuromedin L; Substance K; Neuropeptide K; NPK; TAC1; NKA, NKNA, TAC2; Tachykinin 1; tachykinin 2;
Accession No.	Uniprot:P20366
Calculated MW	kDa
Concentration	1.0 mg/ml
Formulation	Rabbit IgG in 10mM phosphate buffered saline , pH 7.4, 150mM sodium chloride, 0.05% BSA, 0.02% sodium azide and 50% glycerol.
Storage	Store at +4°C for short term. Store at -20°C for long term. Avoid freeze/thaw cycle.

Application Details

IHC 1:100-1:200

Images



Immunohistochemical analysis of paraffin-embedded human brain, using Substance P Antibody.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.