

PF4 Rabbit mAb

Catalog No: #56223

Package Size: #56223-1 50ul #56223-2 100ul


Description

Product Name	PF4 Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB IHC IP
Species Reactivity	Human Rat
Specificity	PF4 Antibody detects endogenous levels of total PF4
Immunogen Description	A synthesized peptide derived from human PF4
Conjugates	Unconjugated
Other Names	C-X-C motif chemokine 4; CXCL4; Iroplact; OncostatinA; PF4; Platelet factor 4; SCYB4; short form; Small inducible cytokine subfamily member 4;
Accession No.	Uniprot:P02776
Calculated MW	11kDa
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at +4ℳ C short term. Store at -20ℳ C long term. Avoid freeze / thaw cycle.

Application Details

WB:1:500~1:2000IHC:1:50~1:200IP:1:50

Images



Western blot analysis of PF4 expression in human spleen lysate.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.