

Langerin (CD207) Rabbit mAb

Catalog No: #56269

Package Size: #56269-1 50ul #56269-2 100ul


Description

Product Name	Langerin (CD207) Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	IHC
Species Reactivity	Human
Specificity	Langerin (CD207) Antibody detects endogenous levels of total Langerin (CD207)
Immunogen Description	A synthesized peptide derived from human Langerin (CD207)
Conjugates	Unconjugated
Other Names	CD207; CD207 antigen; CD207 antigen, langerin; CD207 molecule; CLEC4K; Langerhans cell specific c type lectin; Langerin;
Accession No.	Uniprot:Q9UJ71
Calculated MW	kDa
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at -20°C for long term preservation (recommended). Store at 4°C for short term use.

Application Details

IHC:1:50~1:200

Images



Immunohistochemical analysis of paraffin-embedded human skin, using Langerin (CD207) Antibody.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.