

Bub1 Rabbit mAb

Catalog No: #56369

Package Size: #56369-1 50ul #56369-2 100ul


Description

Product Name	Bub1 Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB IHC
Species Reactivity	Human Mouse Rat
Specificity	Bub1 Antibody detects endogenous levels of total Bub1
Immunogen Description	A synthesized peptide derived from human Bub1
Conjugates	Unconjugated
Other Names	Bub1; BUB1A; BUB1L; hBUB1;
Accession No.	Uniprot:O43683
Calculated MW	122kDa
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at +4ℳ C short term. Store at -20ℳ C long term. Avoid freeze / thaw cycle.

Application Details

WB:1:500~1:1000IHC:1:50~1:200

Images



Western blot analysis of Bub1 expression in K562 cell lysate.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.