

## TCL1 Rabbit mAb

Catalog No: #56383

Package Size: #56383-1 50ul #56383-2 100ul


## Description

|                       |                                                                                                    |
|-----------------------|----------------------------------------------------------------------------------------------------|
| Product Name          | TCL1 Rabbit mAb                                                                                    |
| Host Species          | Rabbit                                                                                             |
| Clonality             | Monoclonal                                                                                         |
| Isotype               | Rabbit IgG                                                                                         |
| Purification          | Affinity-chromatography                                                                            |
| Applications          | WB IHC ICC/IF FC                                                                                   |
| Species Reactivity    | Human                                                                                              |
| Specificity           | TCL1 Antibody detects endogenous levels of total TCL1                                              |
| Immunogen Description | A synthesized peptide derived from human TCL1                                                      |
| Conjugates            | Unconjugated                                                                                       |
| Other Names           | Anti TCL1A; P14 TCL1; T cell leukemia 1; TCL1; TCL1 PEN; Tcl1a;                                    |
| Accession No.         | Uniprot:P56279                                                                                     |
| Calculated MW         | 13kDa                                                                                              |
| Formulation           | Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol. |
| Storage               | Store at +4ℳ C short term. Store at -20ℳ C long term. Avoid freeze / thaw cycle.                   |

## Application Details

WB:1:500~1:2000IHC:1:50~1:200ICC/IF:1:50~1:200FC:1:100

## Images



Western blot analysis of TCL1 expression in Ramos cell lysate.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.