

UCP1+UCP3 Rabbit mAb

Catalog No: #56456

Package Size: #56456-1 50ul #56456-2 100ul


Description

Product Name	UCP1+UCP3 Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB IP
Species Reactivity	Mouse Rat
Specificity	UCP1+UCP3 Antibody detects endogenous levels of total UCP1+UCP3
Immunogen Description	A synthesized peptide derived from human UCP1+UCP3
Conjugates	Unconjugated
Other Names	UCP1; UCP3;
Accession No.	Uniprot:P12242/P56501
Calculated MW	34kDa
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at +4ℳ C short term. Store at -20ℳ C long term. Avoid freeze / thaw cycle.

Application Details

WB:1:500~1:1000IP:1:50

Images



Western blot analysis of UCP1+UCP3 expression in mouse skeletal muscle cell lysate.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.