

ATP2A1/SERCA1 Rabbit mAb

Catalog No: #56484

Package Size: #56484-1 50ul #56484-2 100ul


Description

Product Name	ATP2A1/SERCA1 Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB IHC IP
Species Reactivity	Human
Specificity	ATP2A1/SERCA1 Antibody detects endogenous levels of total ATP2A1/SERCA1
Immunogen Description	A synthesized peptide derived from human ATP2A1/SERCA1
Conjugates	Unconjugated
Other Names	ATP2A; ATP2A1; SERCA1;
Accession No.	Uniprot:O14983
Calculated MW	100kDa
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at +4ℳ C short term. Store at -20ℳ C long term. Avoid freeze / thaw cycle.

Application Details

WB:1:1000~1:5000IHC:1:50~1:200IP:1:50

Images



Western blot analysis of ATP2A1/SERCA1 expression in human fetal muscle lysate.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.