

Syntaxin Rabbit mAb

Catalog No: #56510

Package Size: #56510-1 50ul #56510-2 100ul


Description

Product Name	Syntaxin Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB IP
Species Reactivity	Human Mouse Rat
Specificity	Syntaxin Antibody detects endogenous levels of total Syntaxin
Immunogen Description	A synthesized peptide derived from human Syntaxin
Conjugates	Unconjugated
Other Names	EPIM; Epimorphin; EPM; HPC 1; P35 1; STX1A; Stx1b; STX1B1; STX1B2; STX2; STX2A; STX2B; STX2C; STX3; STX3A;
Accession No.	Uniprot:P61266
Calculated MW	33kDa
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at +4ℳ C short term. Store at -20ℳ C long term. Avoid freeze / thaw cycle.

Application Details

WB:1:1000~1:5000IP:1:50

Images



Western blot analysis of Syntaxin expression in Human fetal brain lysate.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.